

TEKNA ADVANCE SWING



Seat width	36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Front from 40 a 58 cm (adjustable) Rear from 38 a 46 cm (adjustable)
Footplate distance	From 35 to 54 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	95° - 110° - 120°
Total length	85 cm approx.
Total width	Seat width + 17 cm (with 0° camber)
Encumbrance	28 cm approx. (with 0° camber)
Backrest angle	84° - 90° - 96° (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 8,6 kg approx.
Wheelchair weight	Starting from 11,2 kg approx.
Maximum load	125 kg