

# Triride®

by Gianni Conte

A NEW FREEDOM



## ACCESSORIES & SPARE PARTS CATALOG





**IBS** 2  
INTELLIGENT BRAKING SYSTEM

**ACC** 4  
ADAPTIVE CRUISE CONTROL

**ACCESSORIES** 6

**CUSTOMIZATIONS** 10

**SPARE PARTS** 14


**WHEELCHAIR'S TYRE** 18

# IBS Intelligent Braking System

Made in Triride

Triride presents the new Intelligent Braking System that goes beyond the traditional electrical brake. The de inite solution in order to cancel any mechanical and hydraulic brake.

**IBS is an electronic braking system, programmable and customizable by the exclusive Triride software.**



*IBS makes driving very comfortable and easy, improving the stability and safety also thanks to the new anti-lock system (ABS).*

*IBS allows to set up the braking on user's requirements, different surfaces and different driving styles. And all this without using the mechanical brake.*

- **No more brakes pads wear**
- **Exponential increasing of battery charging during braking**
- **Full customization**



*IBS can be used by brake lever or by the electronic brake button that emphasizes the IBS function excluding totally the mechanical brake.*

**YOU MUST HAVE IT ON YOUR  
TRIRIDE!**



# ACC Adaptive Cruise Control

Made in Triride

*Your new technological assistant for all types of routes.  
The advanced electronics that allows the automatic maintenance of the  
desired set speed along all kind of paths in total tranquility and efficiency.*

*The New Adaptive Cruise Control always maintains the set speed,  
even downhill!*



AUTOMATIC MAINTENANCE OF THE DESIRED SET SPEED

## UPHILL

*Automatically  
acting on the  
motor  
acceleration,  
adap-ting itself  
to the slopes of  
the terrain.*

## FLAT

*Acting  
constantly on  
the motor  
acceleration,  
using the  
minor energy  
possible.*

## DOWNHILL

*Maintaining the  
set speed with  
automatic  
activation of  
engine brake.*

**Self-charging  
battery.**

**ACC - A TRIRIDE REVOLUTION FOR DRIVING THE ELECTRICAL DEVICES** *Triride development and research are constantly evolving to implement security and accessibility of your mobility device. The revolutionary Adaptive Cruise Control system is an example. Great autonomy in the gradients for an even more comfortable ride, also reducing the energy consumption thanks to the battery recharging downhill.*

- **Authentic driving assistance both uphill and downhill facilitating the guidance on these paths to the less experienced.**
- **Relax your hand from the accelerator and therefore face the paths in total comfort and carefree.**
- **Operated simply by pressing the dedicated button while driving. Simply deactivation using the brake or the accelerator.**
- **Auto-recharge battery on downhill**



**YOU MUST HAVE IT ON YOUR TRIRIDE!**





BS-O-115

Triride Carrying Bag with Trolley (mod. Foldable)



BS-O-118 / BS-O-119 / BS-O-213

Triride Carrying bag - Size M/L/XL



B-O-211

Bike basket



BS-O-214

Waterproof legs cover



BS-R-346 / BS-R-347

Brackets Neoprene / Neoprene with logo.



BS-O-229

Rear-view mirror



BS-O-172

Neoprene Battery Cover



BS-O-405

Standard front light



BS-O-232

1200 Lumen front light



BS-O-198

USB front light



BS-O-180

USB rear light



BS-O-177

Long fender for 14" wheel



**BS-O-216**

Composite fender



- BS-O-245 black
- BS-O-166 lemon green
- BS-O-167 royal blue
- BS-O-168 pale blue
- BS-O-169 orange

**Parafango in PP**

Standard coloured fender



**BS-O-215**

Full Carbon fender



- BS-O-301 black
- BS-O-127 orange
- BS-O-128 royal blue
- BS-O-129 red
- BS-O-130 lemon green

**Vite a manopola**

Screw knob 8mm x 16mm (couple)



**BS-O-251**

Small Truck trailer



- BS-O-241 Ø 25
- BS-O-116 Ø 30
- BS-O-112 Ø 32
- BS-O-117 Ø 35
- BS-O-111 Ø 40

**BS-O-241**

Extra binding for Truck trailer



BS-O-171

NIK - Parking and cleaning stand



**BS-O-246**

Straight handlebar 42 cm - matt black



**BS-O-208**

4 folds handlebar - matt black



**BS-O-107**

4 folds handlebar 46 cm - matt black



**BS-O-193**

Folding handlebar



**BS-O-107+BS-O-308**

Handlebar with adjustable bar-end  
for tetraplegic users



**BS-O-175**

Mono Command for Tetraplegic Users  
(brake / accelerator)



BS-O-222 / BS-O-220

Adjustable handlebar connection 110/130 mm



BS-O-185 / BS-O-211

Adjustable handlebar connection 160/190 mm



BS-R-299

Thumb accelerator (ambidextrous)



S-R-328 (dx) - S-R-329 (sx)

Thumb accelerator (right - left)



BS-R-364

2 functions control (Base, Foldable, L'8)



BS-R-365

3 functions control (Special)



**BS-O-186**

Multifunctional button



**BS-R-368**

Cruise control button



**BS-R-366**

Electronical brake button



**BS-O-203(sx) / BS-O-273 (dx)**

Brake with locking lever



**BS-O-274(sx) / BS-O-275 (dx)**

Hydraulic brake



**BS-O-235/236/237**

Couple of brackets 29cm / 39cm / customizable



BS-O-272

Curved brackets for some wheelchair models



BS-O-421 (dx) / BS-O-422 (sx)

Up and down brackets (right - left)



BS-O-259 + BS-O-201

Stand with tilting frame





**BS-C-179 / BS-C-180**

Bindings for Mono Wheelchair 70/90



**BS-O-262**

Tilting frame + brackets for Mono Wheelchair



**Consultare il listino per i codici - Check price list for article codes**

Round, oval, elliptical bindings. Available for every kind of wheelchair.



**BS-C-114**

Ottobock Avantgarde CV wheelchair bindings



**BS-C-115**

Offcar Vega wheelchair bindings



**BS-C-118**

Removable footrests wheelchairs "L" bindings (for some wheelchair model only)



**BS-C-222**

Progeo Carbomax wheelchair bindings



**BS-C-148**

Ø25 Basic light bindings for wheelchairs with retractable footrests. (for some wheelchair model only)



**BS-C-085**

Progeo NOIR 2.0 wheelchair bindings



**BS-C-084**

Progeo DUKE wheelchair bindings



**BS-O-218**

Double brake with 2 brake calipers



**BS-R-329 / BS-R-348**

Mechanical brake caliper for Ø180 / Ø160 brake disk



**BS-R-321**

Couple brake pads (Shimano)



**BS-R-358**

Mechanical Couple brake pads (Artek)



**BS-R-357**

Hydraulic couple brake pads DS-17



**BS-C-508**

Tilting frame with stand binding



**BS-O-238**

Spring cotter pin



**BS-R-330**

Screw handle 8mm x 16mm (couple)



**BS-R-331**

Screw handle 8mm x 20mm



**BS-C-186**

Display LCD 01



**BS-C-076**

Display LCD 05



**BS-C-077**

Display LCD 07 with USB port



**BS-R-344**

Cylindrical battery 36V x 7Ah (mod. Base, Foldable, L'8)



**BS-R-313**

Vertical extraction 36V battery (mod. Base)



**BS-R-314 / BS-R-315**

Front extraction 36V/48V battery (Special models)



**BS-C-510 / BS-C-511**

Side extraction 36V/48V battery (mod. MadMax)



**BS-R-310 / BS-R-311 / BS-R-322**

36V / 48V Battery Charger



**BS-R-300(12") / BS-R-303(14")**

Run flat inner tube



**BS-R-339 (12") / BS-R-304 (14") / BS-R-340 (16")**

Inner tube



**BS-R-297 (12") - BS-R-302 (14") - BS-R-298 (16")**

Cross Tyre



**BS-R-001 (12")**

3.0 - 8 Tyre (for Special Light Model)



**BS-R-353 (12") - BS-R-301 (14") - BS-R-334 (16")**

City Tyre



CA-R-100 / CA-R-114

Inner tube 24" x 1,75" / 24" x 1,75"



CA-R-126

Complete Fat wheel 20" x 4"



CA-R-127

Complete Country wheel 24" x 1,75"



CA-R-128

Schwalbe Tyre Marathon Plus 24"x1,75"



CA-R-129

Complete Dynamic wheel 24" x 1,75" (18 rays)



CA-R-130

Schwalbe Marathon 24" x 1 3/8" for wheel 24" x 1"



CA-R-131

Schwalbe Tyre Marathon 24" x 1,75"



CA-R-132

Schwalbe Tyre Marathon 24" x 1"



CA-R-133

Rightron 24" x 1"



CA-R-134

Schwalbe Tyre Road Cruiser 24" x 1,75"



## TRIRIDE IN EUROPE AND IN THE WORLD

Triride srl distributes its products around the World through a commercial network of specialized companies. You can contact our sales representatives in order to test the product and have all the information in the following countries:

- AUSTRIA
  - BALTIC STATES
  - BELGIUM
  - BULGARIA
  - DENMARK
  - FRANCE
  - GERMANY
  - GREECE
  - HUNGARY
  - LUSSEMBOURG
  - NETHERLANDS
  - SPAIN
  - SWEDEN
  - SLOVENIA
  - SWITZERLAND
  - UNITED KINGDOM
- 
- AUSTRALIA
  - CANADA
  - NEW ZEALAND
  - SOUTH AFRICA
  - UNITED ARAB EMIRATES

**FIND OUT YOUR NEAREST RESELLER**  
**RECARE - INFO@RECARE.CO.UK**

*Triride*<sup>®</sup>  
by Gianni Conte

Copyright © 2018 TRIRIDE srl. - Via Massimo D'Antona, 8 – 63812 Montegranaro (FM) - ITALY  
Tel: +39 0733 801405 - Tel/Fax: +39 0733 896964 - Mob. +39 327 6266267 - Email: info@recare.co.uk

[WWW.TRIRIDEITALIA.COM](http://WWW.TRIRIDEITALIA.COM)    Triride