

A NEW FREEDOM







Medical Devices Register of the Health Ministry. TUV Certification No. 15-PEL-0010022-TIC





36 V



540 W



7,5 Kg⁽¹⁾



≥ 40 Km⁽²⁾

TECHNICAL

FEATURES	
Max Power 540W limited & controlled	LCD Display
Electronics by Triride	36V Lithium-ion battery
Reverse	Powerful disc brake
Cruise Control	Tire selection: City or Cross
Assisted Load	Stainless steel frame
Electronic Brake	Speed compliant with legal standards
Weight: 7,5 Kg	5 Speed levels
12" driving wheel	Customized frame color

¹ Weight refers to device excluding battery and arms for wheelchair connection.
² Autonomy refers to device with EcoDrive installed, 36 V 10 Ah battery on a flat route at medium speed of 6 km/h (level 1) with user weight 70 kg. The autonomy depends on user weight, road conditions, battery model, speed/driving mode and on the general state of the wheelchair and tires. Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.