

Medical Devices Register of the Health Ministry. TUV Certification No. 15-PEL-0010022-TIC







36 V / 48 V



750 W / 1000 W



9 Kg⁽¹⁾

Air Chamber 3,0-8 - valve 90



≥ 50 Km⁽²⁾

TECHNICAL FEATURES Max Power 36V 750W / 48V 1000W Electronics by Triride Reverse + Assisted Load Cruise Control Electronic Brake + Powerful disc brake 12" or 14" aluminium driving wheel LCD Display with 5 level speed Tire selection: City or Cross Stainless steel frame Stainless steel frame Customized frame color ALLROAD VERSION 12" - 3.00 - 8 Cross Tyre

36V or 48V Lithium-ion battery

¹Weight refers to device excluding battery and arms for wheelchair connection.

² Autonomy refers to device with EcoDrive installed, 48 V battery on a flat route at medium speed of 6 km/h (level 2) with user weight 70 kg. The autonomy depends on user weight, road conditions, battery model, speed/driving mode and on the general state of the wheelchair and tires.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.