



## TRIBIKE HYBRID PLUS

£4,500.00

\*\*\*

ORDER

☐

QUOTE

☐

ORDER REFERENCE

DATE

Please complete and send this form by e-mail to: [orders@triride.uk](mailto:orders@triride.uk)

### TECHNICAL FEATURES

- |  |  |
|--|--|
| 1- Motor Power 250W  | 8- Shimano 7-speed gearbox   |
| 2- Electronics by Triride  | 9- Disc brake + V-BRAKE brake  |
| 3- Bluetooth KIT: accelerator, brake, reverse, motor brake, battery recharge, IBS, ACC | 10- Speed with engine support max 25   |
| 4- Battery 20-19 with USB: 36V or 48V  | 11-Weight: 11 Kg approx.(without battery)  |
| 5- Pedal Assist PAS  | 12-Adjustable pedaling position - N.B.: after each adjustment check the chain length |
| 6- Display LCD with 5 assistance levels  |  |
| 7- 20" aluminium rim - double chamber  |  |

#### USER'S PERSONAL DATA\*

Name

Surname

Phone

E-mail

Weight

Age

Pathology

#### DELIVERY ADDRESS

Delivery to

Street

Town

Postcode

#### WHEELCHAIR FEATURES

Brand

Model

Front tubes diameter  mm

Seat width  cm

Seat depth (b)  cm

Toe to front tube distance (a)\*  cm


Ground/knee height (d)  cm

Length of stretched arms from the back(c)  cm

Footrest

Fixed	<input type="checkbox"/>	Internal	<input type="checkbox"/>
Removable	<input type="checkbox"/>	External	<input type="checkbox"/>

\* For wheelchairs with removable footrests, measure as per "a" from toe to front tube.



CONFIGURATION				
Wheel	16" wheel		<input type="checkbox"/>	standard
	20" wheel		<input type="checkbox"/>	standard
Battery	BS-R-394	Battery 36V 8.8 USB (316Wh)	<input type="checkbox"/>	standard
	BS-R-393-S	Battery 36V 14 USB (504 Wh)	<input type="checkbox"/>	£270.00
	BS-R-395-S	Battery 48V 10.5 USB (504 Wh)	<input type="checkbox"/>	£270.00
Crancks	TR-B-021	140mm (recommended for 16" wheels )	<input type="checkbox"/>	standard
	TR-B-004	155mm (recommended for 20" wheels)	<input type="checkbox"/>	standard
Crown FSA	TR-B-028	42 cogs (recommended Tribike Hybrid)	<input type="checkbox"/>	standard
	TR-B-027	38 cogs (recommended Tribike Manual)	<input type="checkbox"/>	standard
Handlebar Stem	BS-O-221	Fixed 60mm	<input type="checkbox"/>	standard
	TR-B-090	Fixed 80mm 40 degrees	<input type="checkbox"/>	£35.00
	TR-B-091	Fixed 100mm 40 degrees	<input type="checkbox"/>	£35.00
	BS-C-441	Fixed 120mm 40 degrees	<input type="checkbox"/>	£35.00
Tyre	CA-R-120	City	<input type="checkbox"/>	standard
	CA-R-144	Cross	<input type="checkbox"/>	standard
	CA-R-010	Marathon Schwalbe	<input type="checkbox"/>	£14.00
Colour	BS-O-217-N	Black Ral 9005	<input type="checkbox"/>	standard
	BS-O-217-R	Red Ral 3020	<input type="checkbox"/>	standard
	BS-O-217-W	White Ral 9003	<input type="checkbox"/>	£132.00
	BS-O-217-G	Yellow Ral 1023	<input type="checkbox"/>	£132.00
	BS-O-217-A	Orange Ral 2004	<input type="checkbox"/>	£132.00
	BS-O-217-B	Blue Ral 5015	<input type="checkbox"/>	£132.00
	BS-O-217	Customized RAL RAL no. _____	<input type="checkbox"/>	£177.00
OPTIONAL				
T-Stand with wheels		BS-O-201	<input type="checkbox"/>	£109.00
Leg cover		BS-O-214	<input type="checkbox"/>	£15.00
Laser engraving (customized serigraphy on the frame)		BS-O-174	<input type="checkbox"/>	£68.00
Light - Front 1200 Lumen		BS-O-232	<input type="checkbox"/>	£100.00
Light - Front+USB		BS-O-198	<input type="checkbox"/>	£27.00
Light - Rear+USB		BS-O-180	<input type="checkbox"/>	£32.00
***Prices are intended for the public and before tax				
<b>N. B.:</b> The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product				