

Triride®

by Gianni Conte

A NEW FREEDOM



mod. **TRIBIKE E PLUS Hybrid**



Medical Devices Register of the Health Ministry. TUV Certification No. 15-PEL-0010022-TIC



20"



36 V



250 W



12,5 Kg⁽¹⁾



up to 100 Km⁽²⁾

TECHNICAL FEATURES

EN

36V Lithium-ion battery

Fully adjustable pedal position

Customizable Ergal pedals

PAS assisted riding sensor

Programmable and customizable assistance

Electronics by Triride

USB LCD Display with 5 levels of assistance

Double brake (Disk brake e V-brake)

7-speed Shimano Gear

Tyre: City or Cross

Steel customizable color frame

BLUETOOTH SYSTEM:

IBS - ACC - Reverse

Advance Battery Regen

Start & Walk Assist

¹ Weight refers to device excluding battery and arms for wheelchair connection.

² Autonomy refers to 36V battery on a flat route at medium speed of 6 km/h (level 1) with user weight 70 kg. The autonomy depends on user weight, road conditions, battery model, speed/driving mode and on the general state of the wheelchair and tires.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.