

Outdoor front end Overcoming the barriers of everyday life more easily



Outdoor front end for the outside area



Thresholds and other
Overcomes obstacles more easily
Easy to assemble with just three easy steps



Unique people, unique solutions.

Outdoor front end

Available for the following models:

- Mio (Design 2018)
- Mio Carbon (Design 2018)
- Mio Move
- Vector
- Vector DLS
- Jump alpha
- Jump beta (Sport- and DLS-Version)
- Tilty Vario
- Loop^{SORG}
- Loop^{SORG} RS



Overview

- can easily be mounted and demounted
- all of our wheelchairs can also be equipped after order with an outdoor front end
- splash guard optionally available
- can be extended in width and length on growable wheelchairs, without purchasing parts
- the assembly on makes from other manufactures is on demand, in individual cases possible

Mastering everyday life better

Supposedly easy everyday tasks often pose a real challenge for wheelchair drivers: uneven ground, gravel paths, tram rails, potholes, forest tracks, and the like. Even trying to get to the kindergarden may turn into an unpleasant adventure because the small casters are not suitable for the cobblestone pavements you find in many pedestrian areas.

We have created a simple and clever solution, opening up whole new worlds: the Outdoor Front End. It enables the wheelchair user to navigate on their everyday paths a lot more easily, and in many cases even entirely independently. Furthermore, the Outdoor Front End offers an ideal pushing aid with which both the wheelchair user and their attendant may overcome their daily obstacles with significantly less effort.



Equipment examples

The supplied jack up aid facilitates or enables the independent mounting and dismounting of the outdoor front end.



jack up aid



Splash guard