





## Veloce - REAR Seat-to-floor Height


| WHEEL SIZES                              | Height available |                    | Specify desired height:<br>_____ |
|------------------------------------------|------------------|--------------------|----------------------------------|
| 22" (501) wheel <input type="checkbox"/> | 34.3cm / 13.5"   | to 46.4cm / 18.25" |                                  |
| 24" (540) wheel <input type="checkbox"/> | 36.8cm / 14.5"   | to 48.9cm / 19.25" |                                  |
| 25" (559) wheel <input type="checkbox"/> | 38.1cm / 15"     | to 50.2cm / 19.75" |                                  |
| 26" (590) wheel <input type="checkbox"/> | 39.4cm / 15.5"   | to 51.4cm / 20.25" |                                  |

Configuration used: 2.5 cm (1") seat slope with a center of gravity at 5cm (2")

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1,25 cm (½") of selected Seat to Floor height"

| FRONT FRAME ANGLE |      |                       |                                                                                    |  |  |    |                          |
|-------------------|------|-----------------------|------------------------------------------------------------------------------------|--|--|----|--------------------------|
| 22                | FA70 | 70° Front Frame Angle |  |  |  | NC | <input type="checkbox"/> |
| 23                | FA75 | 75° Front frame angle |                                                                                    |  |  | NC | <input type="checkbox"/> |
| 24                | FA80 | 80° Front Frame Angle |                                                                                    |  |  | NC | <input type="checkbox"/> |
| 25                | FA85 | 85° Front frame angle |                                                                                    |  |  | NC | <input type="checkbox"/> |
| 26                | FA90 | 90° Front Frame Angle |                                                                                    |  |  | NC | <input type="checkbox"/> |

| SEAT TO FOOTREST HEIGHT      |  |                                                                                     |  |  |  |
|------------------------------|--|-------------------------------------------------------------------------------------|--|--|--|
| Set Length to <sup>1</sup> : |  | Length must be at least 4,5 cm (1 ¾") shorter than front seat-to-floor height       |  |  |  |
| 1 - +/- 1.27 cm (1/2")       |  |  |  |  |  |

| Footrest width / Taper selection |     |                                 |       |                                  |       |                                                                                      |  |
|----------------------------------|-----|---------------------------------|-------|----------------------------------|-------|--------------------------------------------------------------------------------------|--|
| Seat Width                       |     | Straight <sup>1</sup>           |       | 1" Taper                         |       |  |  |
| Outside Measurement              |     | 6,4 cm Narrower than Seat Width |       | 11,4 cm Narrower than Seat Width |       |                                                                                      |  |
| 35,56                            | 14" | 29,21                           | 11,5" | 24,13                            | 9,5"  |                                                                                      |  |
| 38,1                             | 15" | 31,75                           | 12,5" | 26,67                            | 10,5" |                                                                                      |  |
| 40,64                            | 16" | 34,29                           | 13,5" | 29,21                            | 11,5" |                                                                                      |  |
| 43,18                            | 17" | 36,83                           | 14,5" | 31,75                            | 12,5" |                                                                                      |  |
| 45,72                            | 18" | 39,37                           | 15,5" | 34,29                            | 13,5" |                                                                                      |  |
| 48,26                            | 19" | 41,91                           | 16,5" | 36,83                            | 14,5" |                                                                                      |  |
| 50,8                             | 20" | 44,45                           | 17,5" | 39,37                            | 15,5" |                                                                                      |  |

1- May interfere with auto-folding and flip-up footplates

| FOOTPLATES            |      |                                                                                           |  |  |      |                          |
|-----------------------|------|-------------------------------------------------------------------------------------------|--|--|------|--------------------------|
| 27                    | FP06 | One piece adjustable angle flip-up ( <i>Aluminum</i> )                                    |  |  | NC   | <input type="checkbox"/> |
| 28                    | FP07 | One piece adjustable angle flip-up ( <i>Carbon fiber with rubber inserts</i> )            |  |  | £365 | <input type="checkbox"/> |
| 29                    | FP28 | One piece auto-folding adjustable angle ( <i>Aluminum</i> )                               |  |  | £240 | <input type="checkbox"/> |
| 30                    | FP29 | One piece auto-folding adjustable angle ( <i>Carbon fiber with rubber inserts</i> )       |  |  | £295 | <input type="checkbox"/> |
| 31                    | FP02 | Newton adjustable angle footplate ( <i>Composites - Two piece</i> )                       |  |  | £175 | <input type="checkbox"/> |
| 32                    | FP03 | Newton adjustable angle ( <i>Aluminum - Two piece</i> )                                   |  |  | £192 | <input type="checkbox"/> |
| 33                    | FP04 | Newton adjustable angle footplate ( <i>Carbon fiber with rubber inserts - Two piece</i> ) |  |  | £345 | <input type="checkbox"/> |
| FOOTPLATE ACCESSORIES |      |                                                                                           |  |  |      |                          |
| 34                    | FP08 | Length adjustable heel loop ( <i>not available with one piece footplate</i> )             |  |  | £35  | <input type="checkbox"/> |
| 35                    | FP09 | Calf strap                                                                                |  |  | £30  | <input type="checkbox"/> |
| 36                    | FP30 | Bodypoint Aeromesh Padded calf strap                                                      |  |  | £95  | <input type="checkbox"/> |

| Backrest                                             |      |                                                                                                                                                                                                                                    |                                                      |  |      |                          |
|------------------------------------------------------|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------|--|------|--------------------------|
| FIXED HEIGHT, NO PUSH HANDLE OPTIONS                 |      |                                                                                                                                                                                                                                    |                                                      |  |      |                          |
| 37                                                   | BC03 | Ultralight carbon fiber back cane <sup>1</sup><br>12 (9") <input type="checkbox"/> 22 (12") <input type="checkbox"/> 38 (15") <input type="checkbox"/> 45,7 (18") STD <input type="checkbox"/> 53,3 (21") <input type="checkbox"/> |                                                      |  | NC   | <input type="checkbox"/> |
| 38                                                   | BC04 | 8° bend carbon fiber back cane <sup>1</sup><br>34,3 (13 ½) <input type="checkbox"/> 38,1 (15") <input type="checkbox"/> 45,7 (18") STD <input type="checkbox"/>                                                                    |                                                      |  | £40  | <input type="checkbox"/> |
| HEIGHT ADJUSTABLE WITH PUSH HANDLES                  |      |                                                                                                                                                                                                                                    |                                                      |  |      |                          |
| 39                                                   | BC06 | Straight carbon fiber back cane<br>22 - 30,4 (9-12") <input type="checkbox"/> 30,4 - 38 (12-15") <input type="checkbox"/> 38 - 45,7 (15-18") <input type="checkbox"/> 45,7 - 53,3 (18-21") STD <input type="checkbox"/>            |                                                      |  | NC   | <input type="checkbox"/> |
| 40                                                   | BC08 | 8° bend Carbon fiber back cane<br>34,3 - 38 (13 ½-15") <input type="checkbox"/> 38,1 - 45,7 (15-18") <input type="checkbox"/> 45,7 - 53,3 (18-21") STD <input type="checkbox"/>                                                    |                                                      |  | £40  | <input type="checkbox"/> |
| HEIGHT AND ANGLE ADJUSTABLE WITH PUSH HANDLE OPTIONS |      |                                                                                                                                                                                                                                    |                                                      |  |      |                          |
| 41                                                   | BC09 | Straight carbon fiber back cane<br>30,4 - 38 (12-15") <input type="checkbox"/> 38,1 - 45,7 (15-18") <input type="checkbox"/> 45,7 - 53,3 (18-21") STD <input type="checkbox"/>                                                     |                                                      |  | £160 | <input type="checkbox"/> |
| 42                                                   | BC10 | 8° bend carbon fiber back cane<br>34,3 - 38 (13 ½-15") <input type="checkbox"/> 38,1 - 45,7 (15-18") <input type="checkbox"/> 45,7 - 53,3 (18-21") STD <input type="checkbox"/>                                                    |                                                      |  | £210 | <input type="checkbox"/> |
| 43                                                   |      | Select back cane angle: 85° <sup>1</sup> <input type="checkbox"/> 90° <input type="checkbox"/> 95° STD <input type="checkbox"/> 100° <input type="checkbox"/> 105° <input type="checkbox"/> 110° <input type="checkbox"/>          |                                                      |  |      |                          |
| Set height to: <sup>2</sup>                          |      |                                                                                                                                                                                                                                    | (measured from seat sling to top of back upholstery) |  |      |                          |
| PUSH HANDLE OPTIONS                                  |      |                                                                                                                                                                                                                                    |                                                      |  |      |                          |
| 44                                                   | PH01 | Push handles                                                                                                                                                                                                                       |                                                      |  | NC   | <input type="checkbox"/> |
| 45                                                   | PH05 | Fold down push handles <sup>3, 5</sup>                                                                                                                                                                                             |                                                      |  | £100 | <input type="checkbox"/> |
| 46                                                   | PH06 | Clamp on height adjustable push handle <sup>3, 4, 5</sup> 22,9 -50,8 cm (9-20") STD <input type="checkbox"/> 22,9 - 58.4 cm (9-23") <input type="checkbox"/>                                                                       |                                                      |  | £250 | <input type="checkbox"/> |
| 47                                                   | PH02 | Omit push handles                                                                                                                                                                                                                  |                                                      |  | NC   | <input type="checkbox"/> |

1- Not available with a 8 degree bend cane    2 - If not indicated, default height setting is at the lowest position    3 - 120 kg/265 pounds weight limit

4 - Not compatible with swing away armrests, 8° bend back canes

5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

| BACK UPHOLSTERY |      |                                                |  |  |      |                          |
|-----------------|------|------------------------------------------------|--|--|------|--------------------------|
| 48              | BH07 | Slip-on Active-Mesh back upholstery            |  |  | NC   | <input type="checkbox"/> |
| 49              | BH08 | Adjustable tension Active-Mesh back upholstery |  |  | £140 | <input type="checkbox"/> |
| 50              | BH99 | Omit Back upholstery                           |  |  | NC   | <input type="checkbox"/> |

| SEAT CUSHION |      |                     |  |  |      |                          |
|--------------|------|---------------------|--|--|------|--------------------------|
| 51           | SC01 | 5,1 cm (2") Cushion |  |  | £140 | <input type="checkbox"/> |
| 52           | SC02 | 7,6 cm (3") Cushion |  |  | £155 | <input type="checkbox"/> |
| 53           | SC99 | Omit seat cushion   |  |  | STD  | <input type="checkbox"/> |

| Rear wheels                                                                                                  |      |                                                            |                                         |                               |                                |                               |                              |                              |                               |                              |                              |    |                          |
|--------------------------------------------------------------------------------------------------------------|------|------------------------------------------------------------|-----------------------------------------|-------------------------------|--------------------------------|-------------------------------|------------------------------|------------------------------|-------------------------------|------------------------------|------------------------------|----|--------------------------|
| ONE ARM DRIVE <sup>1</sup> Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims |      |                                                            |                                         |                               |                                |                               |                              |                              |                               |                              |                              |    |                          |
| 54                                                                                                           | HR22 | One arm drive <sup>2, 3</sup>                              | 24"                                     | left <input type="checkbox"/> | right <input type="checkbox"/> |                               |                              |                              |                               | £1 165                       | <input type="checkbox"/>     |    |                          |
| 55                                                                                                           | HR98 | Omit handrim on opposite side of the One Arm Drive (OAD)   |                                         |                               |                                |                               |                              |                              |                               | NC                           | <input type="checkbox"/>     |    |                          |
| REAR WHEELS                                                                                                  |      |                                                            |                                         |                               |                                |                               |                              |                              |                               |                              |                              |    |                          |
| 56                                                                                                           | RW01 | Mag <sup>4</sup>                                           | 22" <input type="checkbox"/>            | 24" <input type="checkbox"/>  |                                |                               |                              |                              |                               | £40                          | <input type="checkbox"/>     |    |                          |
| 57                                                                                                           | RW02 | Newton One - spoke wheel <sup>5</sup>                      | 22" <input type="checkbox"/>            | 24" <input type="checkbox"/>  |                                |                               |                              |                              |                               | NC                           | <input type="checkbox"/>     |    |                          |
| 58                                                                                                           | RW03 | Newton Gravity - ultralight wheel                          | 22" <input type="checkbox"/>            | 24" <input type="checkbox"/>  | 25" <input type="checkbox"/>   | 26" <input type="checkbox"/>  |                              |                              |                               | £250                         | <input type="checkbox"/>     |    |                          |
| 59                                                                                                           | RW04 | Spinerigy Spox                                             | 22" <input type="checkbox"/>            | 24" <input type="checkbox"/>  | 25" <input type="checkbox"/>   | 26" <input type="checkbox"/>  |                              |                              |                               | £590                         | <input type="checkbox"/>     |    |                          |
| 60                                                                                                           |      | Spinerigy Spox spokes colors <sup>6</sup>                  | <b>Blk STD</b> <input type="checkbox"/> | Wht <input type="checkbox"/>  | Blue <input type="checkbox"/>  | Red <input type="checkbox"/>  | Yel <input type="checkbox"/> |                              |                               | NC                           | <input type="checkbox"/>     |    |                          |
| 61                                                                                                           | RW05 | Spinerigy LX                                               | 22" <input type="checkbox"/>            | 24" <input type="checkbox"/>  | 25" <input type="checkbox"/>   | 26" <input type="checkbox"/>  |                              |                              |                               | £665                         | <input type="checkbox"/>     |    |                          |
| 62                                                                                                           |      | Spinerigy LX spokes colors <sup>6</sup>                    | <b>Blk STD</b> <input type="checkbox"/> | Pur <input type="checkbox"/>  | Wht <input type="checkbox"/>   | Blue <input type="checkbox"/> | Red <input type="checkbox"/> | Yel <input type="checkbox"/> | Pink <input type="checkbox"/> | Grn <input type="checkbox"/> | Org <input type="checkbox"/> | NC | <input type="checkbox"/> |
| 63                                                                                                           | RW07 | Spinerigy CLX Carbon Blade                                 |                                         | 24" <input type="checkbox"/>  | 25" <input type="checkbox"/>   |                               |                              |                              |                               | £1 235                       | <input type="checkbox"/>     |    |                          |
| 64                                                                                                           |      | Spinerigy CLX spokes colors <sup>6</sup>                   | <b>Blk STD</b> <input type="checkbox"/> | Pur <input type="checkbox"/>  | Wht <input type="checkbox"/>   | Blue <input type="checkbox"/> | Red <input type="checkbox"/> | Yel <input type="checkbox"/> | Pink <input type="checkbox"/> | Grn <input type="checkbox"/> | Org <input type="checkbox"/> | NC | <input type="checkbox"/> |
| TIRES                                                                                                        |      |                                                            |                                         |                               |                                |                               | Available                    |                              |                               |                              |                              |    |                          |
| 65                                                                                                           | TI01 | Soft Urethane 3,17cm ( 1 1/4")                             |                                         |                               |                                |                               | 22" 24"                      |                              |                               | NC                           | <input type="checkbox"/>     |    |                          |
| 66                                                                                                           | TI02 | Soft Urethane knobby 3,49cm (1 3/8")                       |                                         |                               |                                |                               | 22", 24", 25", 26"           |                              |                               | £80                          | <input type="checkbox"/>     |    |                          |
| 67                                                                                                           | TI04 | Soft Urethane - Shox 2,54cm (1")                           |                                         |                               |                                |                               | 22", 24", 25", 26"           |                              |                               | £96                          | <input type="checkbox"/>     |    |                          |
| 68                                                                                                           | TI05 | Pneumatic 3,49cm (1 3/8")                                  |                                         |                               |                                |                               | 22", 24", 25", 26"           |                              |                               | NC                           | <input type="checkbox"/>     |    |                          |
| 69                                                                                                           | TI08 | Pneumatic with Airless insert 3,49cm (1 3/8") <sup>7</sup> |                                         |                               |                                |                               | 22", 24", 26"                |                              |                               | £66                          | <input type="checkbox"/>     |    |                          |
| 70                                                                                                           | TI06 | Pneumatic HP - 2,54cm (1")                                 |                                         |                               |                                |                               | 22", 24", 25", 26"           |                              |                               | £80                          | <input type="checkbox"/>     |    |                          |
| 71                                                                                                           | TI07 | Pneumatic HP - Schwalbe Marathon Plus 2,54cm (1")          |                                         |                               |                                |                               | 24", 25", 26"                |                              |                               | £90                          | <input type="checkbox"/>     |    |                          |

1- Offered with handrims on each side of the wheelchair
2- Only standard handrim position available
3- Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position
5 - Not compatible with Shox tires

6 - Colour spokes require an additionnal 10 days lead-time
7- Only available with Newton One and Newton Gravity rear wheels

| Handrims |      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |  |  |  |  |     |      |                          |
|----------|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|--|-----|------|--------------------------|
| HANDRIM  |      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |  |  |  |  |     |      |                          |
| 72       | HR01 | Aluminum anodized                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |  |  |  |  |     | STD  | <input type="checkbox"/> |
| 73       | HR61 | Aluminum black hard anodised <sup>1,3</sup>                                                                                                                                                                                                                                                                                                                                                                                                                                                    |  |  |  |  | New | £50  | <input type="checkbox"/> |
| 74       | HR02 | Plastic coated <sup>2</sup>                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |  |  |  |  |     | £80  | <input type="checkbox"/> |
| 75       | HR03 | Q-Grip High friction Neoprene coated <sup>2,3,4</sup>                                                                                                                                                                                                                                                                                                                                                                                                                                          |  |  |  |  |     | £295 | <input type="checkbox"/> |
| 76       | HR07 | Newton Air Grip <sup>3</sup> <b>Black</b> <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                                                                                                                             |  |  |  |  |     | £80  | <input type="checkbox"/> |
| 77       | HR   | Newton Air Grip <sup>3</sup> (08) Sunkissed Orange <input type="checkbox"/> (09) Burgundy <input type="checkbox"/> (10) Sapphire Blue <input type="checkbox"/> (11) Monster Green <input type="checkbox"/> (12) Charcoal <input type="checkbox"/> (13) Fuchsia <input type="checkbox"/> (14) Steel Blue <input type="checkbox"/> (15) Ferrari Red <input type="checkbox"/> (16) Acid Green <input type="checkbox"/> (17) White <input type="checkbox"/> (18) Sky Blue <input type="checkbox"/> |  |  |  |  |     | £200 | <input type="checkbox"/> |
| 78       | HR28 | Surge <sup>2,3,4</sup> (1-1/2" Oval handrim w/ Gripton high friction strip)                                                                                                                                                                                                                                                                                                                                                                                                                    |  |  |  |  |     | £315 | <input type="checkbox"/> |
| 79       | HR05 | Surge LT <sup>2,3,4</sup> (1-1/8" Oval handrim w/ Gripton high friction strip)                                                                                                                                                                                                                                                                                                                                                                                                                 |  |  |  |  |     | £338 | <input type="checkbox"/> |
| 80       | HR04 | Natural Fit Standard Grip <sup>2,3,4</sup>                                                                                                                                                                                                                                                                                                                                                                                                                                                     |  |  |  |  |     | £260 | <input type="checkbox"/> |
| 81       | HR99 | Omit handrim                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |  |  |  |  |     | NC   | <input type="checkbox"/> |

1 - Not available in 26"
2- Only standard handrim position available
3 - Not available with Mag wheels or OAD
4 - Not available in 20"

| Handrims option                                             |       |                           |                                                    |                                               |  |  |  |      |                          |
|-------------------------------------------------------------|-------|---------------------------|----------------------------------------------------|-----------------------------------------------|--|--|--|------|--------------------------|
| Assembly position for aluminum and Newton AirGrip handrims  |       |                           |                                                    |                                               |  |  |  |      |                          |
| 82                                                          | HR29  | Standard Handrim position | Full tab                                           |                                               |  |  |  | STD  | <input type="checkbox"/> |
| 83                                                          | HR30  | Narrow Handrim position   | Short tab 2 positions                              |                                               |  |  |  | £20  | <input type="checkbox"/> |
| 84                                                          | HR31  | Super Narrow Handrim      | Short tab 2 positions STD <input type="checkbox"/> | Short tab 1 position <input type="checkbox"/> |  |  |  | £20  | <input type="checkbox"/> |
| Thumb piece <sup>4</sup> selection for Natural Fit handrims |       |                           |                                                    |                                               |  |  |  |      |                          |
| 85                                                          | HR23  | Standard grip             |                                                    |                                               |  |  |  | £70  | <input type="checkbox"/> |
| 86                                                          | HR27  | Super grip                |                                                    |                                               |  |  |  | £118 | <input type="checkbox"/> |
| 87                                                          | HR099 | Omit thumb piece          |                                                    |                                               |  |  |  | NC   | <input type="checkbox"/> |


4- Compatible with Surge, Surge LT handrims

| AXLE <sup>1</sup> |      |                    |  |  |  |  |  |      |                          |
|-------------------|------|--------------------|--|--|--|--|--|------|--------------------------|
| 88                | AX01 | Quick release axle |  |  |  |  |  | STD  | <input type="checkbox"/> |
| 89                | AX03 | Quad release axle  |  |  |  |  |  | £105 | <input type="checkbox"/> |

1- Not required with one arm drive

| REAR WHEEL CAMBER |      |           |  |  |  |  |  |      |                          |
|-------------------|------|-----------|--|--|--|--|--|------|--------------------------|
| 90                | CM00 | 0° camber |  |  |  |  |  | STD  | <input type="checkbox"/> |
| 91                | CM03 | 3° camber |  |  |  |  |  | NC   | <input type="checkbox"/> |
| 92                | CM06 | 6° camber |  |  |  |  |  | £125 | <input type="checkbox"/> |

# CENTER OF GRAVITY



## Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

X =

| Most stable              |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          | Most responsive          |  |  |  |
|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--|--|--|
| 2,5<br>1"                | 3,2<br>1/4"              | 3,8<br>1/2"              | 4,4<br>3/4"              | 5,1<br>2"                | 5,7<br>1/4"              | 6,3<br>1/2"              | 7<br>3/4"                | 7,6<br>3"                | 8,3<br>1/4"              | 8,9<br>1/2"              | 9,5<br>3/4"              | 10,2<br>4"               |  |  |  |
| <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |  |  |  |

**If not indicated, default setting is at 3,8 cm/1.5"**

Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6,3 cm (2.5").

Possible conflict between a 26" rear wheel and a center of gravity over 6,3 cm (2.5").

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

## WHEEL LOCK

|    |              |                                                                                |  |  |  |      |                          |
|----|--------------|--------------------------------------------------------------------------------|--|--|--|------|--------------------------|
| 93 | WL08         | Newton Wheel lock - push to lock                                               |  |  |  | STD  | <input type="checkbox"/> |
| 94 | WL02         | Pull to lock                                                                   |  |  |  | NC   | <input type="checkbox"/> |
| 95 | WL04         | Newton grade aid push to lock (suggested with pneumatic or high pressure tire) |  |  |  | £125 | <input type="checkbox"/> |
| 96 | WL13<br>WL14 | 6" removable extension handle                                                  |  |  |  | £65  | <input type="checkbox"/> |
| 97 | WL05         | Scissor Lock                                                                   |  |  |  | £120 | <input type="checkbox"/> |
| 98 | WL99         | Omit wheel locks                                                               |  |  |  | NC   | <input type="checkbox"/> |

## Handrims

| FLIP BACK ARMREST <sup>1</sup>  |      |                                              | Height Adjustment     |                              | Removable function     |                                 |                                                                                              |
|---------------------------------|------|----------------------------------------------|-----------------------|------------------------------|------------------------|---------------------------------|----------------------------------------------------------------------------------------------|
| 99                              | AR01 | 10" Desk Length (25,4 cm)                    | 20,3 - 30,5 (8 - 12") | <input type="checkbox"/> STD | 25,4 - 35,6 (10 - 14") | <input type="checkbox"/> Locked | <input type="checkbox"/> Unlocked <input type="checkbox"/> STD £295 <input type="checkbox"/> |
| 100                             | AR02 | 14" Full Length (35,6 cm)                    | 20,3 - 30,5 (8 - 12") | <input type="checkbox"/> STD | 25,4 - 35,6 (10 - 14") | <input type="checkbox"/> Locked | <input type="checkbox"/> Unlocked <input type="checkbox"/> STD £295 <input type="checkbox"/> |
| "T" ARMREST <sup>1</sup>        |      |                                              | Height Adjustment     |                              |                        |                                 |                                                                                              |
| 101                             | AR03 | 10" Desk Length (25,4 cm)                    | 20,3 - 30,5 (8 - 12") | <input type="checkbox"/> STD | 25,4 - 35,6 (10 - 14") |                                 | £270 <input type="checkbox"/>                                                                |
| 102                             | AR04 | 14" Full Length (35,6 cm)                    | 20,3 - 30,5 (8 - 12") | <input type="checkbox"/> STD | 25,4 - 35,6 (10 - 14") |                                 | £270 <input type="checkbox"/>                                                                |
| SWING AWAY ARMREST <sup>2</sup> |      |                                              | Height Adjustment     |                              |                        |                                 |                                                                                              |
| 103                             | AR05 | Tubular arm pad                              | 20,3 - 27,9 (8 - 11") |                              |                        |                                 | £190 <input type="checkbox"/>                                                                |
| 104                             | AR06 | 10 "Desk length arm pad                      | 20,3 - 27,9 (8 - 11") |                              |                        |                                 | £190 <input type="checkbox"/>                                                                |
| GEL ARM PAD <sup>3</sup>        |      |                                              |                       |                              |                        |                                 |                                                                                              |
| 105                             | AR12 | 36cm (14") Full Length Arm Pad - Gel Ovation |                       |                              |                        |                                 | £125 <input type="checkbox"/>                                                                |
| 106                             | AR99 | Omit armrest                                 |                       |                              |                        |                                 | NC <input type="checkbox"/>                                                                  |

1 - May have a 1/2 " adjustment restriction on both limits of the range

2 - Could interfere with some rigid back rests

3- Not available with swing away armrest

## SIDE GUARDS (not available with Flip Back and "T" Armrest)

| 107                | SG01 | Plastic side guards (removable)                     |  |  |  | £165 | <input type="checkbox"/> |
|--------------------|------|-----------------------------------------------------|--|--|--|------|--------------------------|
| 108                | SG02 | Carbon fiber side guards (fixed)                    |  |  |  | £210 | <input type="checkbox"/> |
| 109                | SG03 | Carbon fiber side guards (removable)                |  |  |  | £240 | <input type="checkbox"/> |
| 110                | SG06 | Plastic fender side guards (removable) <sup>1</sup> |  |  |  | £225 | <input type="checkbox"/> |
| 111                | SG04 | Carbon fiber fender side guards (removable)         |  |  |  | £415 | <input type="checkbox"/> |
| 112                | SG99 | Omit side guard                                     |  |  |  | NC   | <input type="checkbox"/> |
| SIDE GUARDS OPTION |      |                                                     |  |  |  |      |                          |
| 113                | SG07 | Side guards clip (pair) <sup>2</sup>                |  |  |  | NC   | <input type="checkbox"/> |

1- Not compatible with 20" and 22" rear wheels

2- Not available with fixed carbon fiber side guard

## POSITIONING BELT

|     |      |                                                                 |                        |                          |                      |                          |      |                          |
|-----|------|-----------------------------------------------------------------|------------------------|--------------------------|----------------------|--------------------------|------|--------------------------|
| 114 | SB01 | Velcro adjustable                                               | Standard 88,9 cm (35") | <input type="checkbox"/> | Long 152,40 cm (60") | <input type="checkbox"/> | £40  | <input type="checkbox"/> |
| 115 | SB02 | Auto buckle                                                     | Standard 88,9 cm (35") | <input type="checkbox"/> | Long 152,40 cm (60") | <input type="checkbox"/> | £70  | <input type="checkbox"/> |
| 116 | SB07 | Hardware attachment for 2 point seatbelt (1 pair)               |                        |                          |                      |                          | NC   | <input type="checkbox"/> |
| 117 | SB03 | Additional attachment hardware for 4 point seatbelt             |                        |                          |                      |                          | £15  | <input type="checkbox"/> |
| 118 | SB04 | BodyPoint Hip Belt - 2 point not padded                         | M (142,24 cm/56")      | <input type="checkbox"/> | L (154,94 cm/61")    | <input type="checkbox"/> | £66  | <input type="checkbox"/> |
| 119 | SB05 | BodyPoint Hip Belt - 2 point padded                             | S (129,54 cm/51")      | <input type="checkbox"/> | M (142,24 cm/56")    | <input type="checkbox"/> | £145 | <input type="checkbox"/> |
| 120 | SB06 | BodyPoint Hip Belt - 4 point padded                             | S (129,54 cm/51")      | <input type="checkbox"/> | M (142,24 cm/56")    | <input type="checkbox"/> | £165 | <input type="checkbox"/> |
| 121 | SB08 | BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps | M (73,66 cm/18")       | <input type="checkbox"/> | L (93,98 cm/24")     | <input type="checkbox"/> | £240 | <input type="checkbox"/> |
| 122 | SB99 | Omit positioning belt                                           |                        |                          |                      |                          | STD  | <input type="checkbox"/> |

## HEADREST <sup>1</sup>

|     |      |                                             |                           |                          |                            |                          |                         |                          |      |                          |
|-----|------|---------------------------------------------|---------------------------|--------------------------|----------------------------|--------------------------|-------------------------|--------------------------|------|--------------------------|
| 123 | AC01 | Profiled Headrest                           | S 8.8 X 15.2 cm (3.5"x6") | <input type="checkbox"/> | M 11.4 X 17.8 cm (4.5"x7") | <input type="checkbox"/> | L 14 X 2.3 cm (5.5"x8") | <input type="checkbox"/> | £155 | <input type="checkbox"/> |
| 124 | AC28 | Headrest mounting hardware with tension bar |                           |                          |                            |                          |                         |                          | £245 | <input type="checkbox"/> |

1 - Only available with push handles

| ACCESSORIES |      |                                                      |                              |                                            |                                |                                           |                                |                                 |                                  |                          |
|-------------|------|------------------------------------------------------|------------------------------|--------------------------------------------|--------------------------------|-------------------------------------------|--------------------------------|---------------------------------|----------------------------------|--------------------------|
| 125         | AC03 | Anti tippers                                         |                              |                                            |                                |                                           |                                | £120                            | <input type="checkbox"/>         |                          |
| 126         | AC29 | Swing-away anti-tippers <sup>1</sup>                 |                              | Left <input type="checkbox"/>              | Right <input type="checkbox"/> | New                                       |                                | £115                            | <input type="checkbox"/>         |                          |
| 127         | AC27 | Tip assist <sup>2</sup>                              |                              |                                            |                                |                                           |                                | £50                             | <input type="checkbox"/>         |                          |
| 128         | AC05 | Transit Tie Down (Unoccupied)                        |                              |                                            |                                |                                           |                                | £270                            | <input type="checkbox"/>         |                          |
| 129         | AC06 | WC-19 Transit Tie Down w/Q'straint belt <sup>3</sup> |                              | Short belt 20-41" <input type="checkbox"/> |                                | Long belt 20-46" <input type="checkbox"/> |                                | £100                            | <input type="checkbox"/>         |                          |
| 130         | AC08 | Spoke guards <sup>4</sup>                            |                              | Available: 22", 24", 26"                   |                                |                                           |                                | £100                            | <input type="checkbox"/>         |                          |
| 131         | AC10 | Seat pouch                                           |                              |                                            |                                |                                           |                                | £60                             | <input type="checkbox"/>         |                          |
| 132         | AC26 | BodyPoint Mobility Bag                               |                              |                                            |                                |                                           |                                | £105                            | <input type="checkbox"/>         |                          |
| 133         | AC23 | Newton Backpack                                      |                              |                                            |                                |                                           |                                | £100                            | <input type="checkbox"/>         |                          |
| 134         | AC12 | Velcro Closing Strap                                 |                              |                                            |                                |                                           |                                | £40                             | <input type="checkbox"/>         |                          |
| 135         | AC13 | Tool kit                                             |                              |                                            |                                |                                           |                                | £40                             | <input type="checkbox"/>         |                          |
| 136         | AC16 | Motion Composites Polo Shirt                         | W/S <input type="checkbox"/> | W/M <input type="checkbox"/>               | W/L <input type="checkbox"/>   |                                           |                                |                                 | £32                              | <input type="checkbox"/> |
|             |      |                                                      | M/S <input type="checkbox"/> | M/M <input type="checkbox"/>               | M/L <input type="checkbox"/>   | M/XL <input type="checkbox"/>             | M/XXL <input type="checkbox"/> | M/XXXL <input type="checkbox"/> | M/XXXXL <input type="checkbox"/> |                          |

| CARBON FIBER FINISH |    |                  |                                           |                                       |                                          |                                                |  |  |    |                          |
|---------------------|----|------------------|-------------------------------------------|---------------------------------------|------------------------------------------|------------------------------------------------|--|--|----|--------------------------|
| 137                 | FC | Highlight color: | (10) Ferrari Red <input type="checkbox"/> | (11) White <input type="checkbox"/>   | (12) Black <input type="checkbox"/>      | (13) Sunkissed Orange <input type="checkbox"/> |  |  | NC | <input type="checkbox"/> |
|                     |    |                  | (14) Sky Blue <input type="checkbox"/>    | (15) Fuchsia <input type="checkbox"/> | (16) Acid Green <input type="checkbox"/> |                                                |  |  |    |                          |
| TOTAL               |    |                  |                                           |                                       |                                          |                                                |  |  |    |                          |

**Notes :**