

Hawk Order Form

Fill in and send to Wolturnus - Ausfüllen und an Wolturnus senden - Lødfylt og send til Wolturnus

Subject to Technical Modifications and Printing Errors / Tekniske ændringer og trykfejl
Forbehold for tekniske ændringer og trykfejl

© Wolturnus A/S OF HAWK Winter 2026 v5.0

User information	
User name	
Address	
Zip-code & city	
Tel. / Telefon	
E-Mail	
Weight	Kg.
Height	cm.
Offer Only	☐
Order	☐
Payment Details	
EAN number	
Your PO number	
Invoice recipient	
Phone/Telefon	
E-Mail	
Contact Person	
Name and address	
Phone	
E-Mail	
Delivery Information	
Consultant	☐
User	☐
Dealer	☐
Other address	☐
Consultant / Date	



Custom-built wheelchair for active use Default Equipment

- Tailor-made rigid-frame wheelchair
- Fully welded or adjustable versions available
- Made of 7020 T6 aluminium
- Adjustable seat and back upholstery
- Impact-resistant anodized frame
- Rear wheels with quick-release axles

£2,995



Orderform is only compatible with
Adobe Acrobat Reader which is available at:

PC & Macbook
<https://get.adobe.com/uk/reader/>

Tablet & iPad
<https://www.adobe.com/acrobat/mobile.html>

☒ ENG 

☐ DE 

☐ DK 

☐ FR 

Remarks / Bemærkninger

REMARKS

This combination/option is available

- This combination/option is **not** available

- ! This information is mandatory and must always be filled in

If no price is mentioned in the left column, the option is included in the price mentioned at the front page. Prices excl VAT.

For options or solutions not mentioned in this orderform contact Wolturnus A/S. See contact information at the bottom of the page.

BEMÆRK

Denne kombination/valg er tilgængeligt

- Denne kombination/valg er **ikke** tilgængeligt

- ! Denne information er nødvendig og skal altid udfyldes


Hvis der ingen pris er angivet i venstre kolonne, er valget inkluderet i prisen angivet på forsiden. Priser ekskl. moms.



For valgmuligheder / løsninger, som ikke er nævnt i dette måleskema kontakt venligst Wolturnus A/S. Se kontakt info nederst på siden.

! Model (mandatory)		
1510 ???? ?09 £2,995	☒	Hawk

! Frame colour (mandatory)		
1510????2109	☐	White Powder coat nr. 21
1510????2509	☐	Black Powder coat nr. 25
1510????2909	☐	Grey Powder coat nr. 29


! Secondary colour (mandatory) / Sekundær farve (skal udfyldes)		
Components in secondary colour: Front fork, foot rest bar, baluster for lift-off sideguards and mounting plate for folding back		
6190000040??	☐	Bronze
619000004004	☐	Blue
619000004005	☐	Red
619000004006	☐	Green
619000004007	☐	Black
619000004008	☐	Sliver

! Seat width and seat depth (mandatory)									
		Seat width (SW)							
Seat depth (SD)		35 cm	37,5 cm	40 cm	42,5 cm	45 cm	47,5 cm	50 cm	
Small	35-42 cm	☐	☐	☐	☐	☐	☐	☐	
Medium	38-45 cm	☐	☐	☐	☐	☐	☐	☐	
Large	42-50 cm	☐	☐	☐	☐	☐	☐	☐	
SD: Requested seat depth:		_____ cm							



! Seat height from front SH-f (mandatory)

		Seat height from front* (SH-f) (cm)												
Front wheel		46 cm	46,5 cm	47 cm	47,5 cm	48 cm	48,5 cm	49 cm	49,5 cm	50 cm	50,5 cm	51 cm	51,5 cm	52 cm
Small frame	4"		-	-		-	-	-	-	-	-	-	-	-
	5"	-	-	-		-	-		-	-	-	-	-	-
	6"	-	-	-	-	-	-	-	-	-	-	-	-	-
Medium frame	4"	-	-		-	-		-	-	-	-	-	-	-
	5"	-	-	-	-	-		-	-		-	-	-	-
	6"	-	-	-	-	-	-	-	-	-	-	-	-	-
Large frame	4"	-	-	-	-		-	-		-	-	-	-	-
	5"	-	-	-	-	-	-	-		-	-		-	-
	6"	-	-	-	-	-	-	-	-	-	-	-	-	



SH-f

Castors

96640100???? £0		4" Fiber plastic front wheel 34 mm width (Std.)	Black ☐ Silver ☐ White ☐ Red ☐ Blue ☐
966401253400 £0		5" Fiber plastic front wheel 34 mm width (nur in schwarz)	
966401503400 £0		6" Fiber plastic front wheel 34 mm width (nur in schwarz) (Nur in 52 cm SF-h)	
860020000000 £45		Anti-haircap for fiber plast wheels	


966401003410 £150		4" Aluminum front wheel 34 mm width (Silver)
966401253410 £150		5" Aluminum front wheel 34 mm width (Silver)
966401004610 £395		4" Carbon Front wheel High gloss 30mm width (Black)
966401254610 £395		5" Carbon Front wheel High gloss 30mm width (Black)

966401000??? £275		4" Aluminum front wheel 35 mm width Froglegs Softwheel	Black ☐ Silver ☐ White ☐ Red ☐ Blue ☐
966401250??? £275		5" Aluminum front wheel 35 mm width Froglegs Softwheel	
966401500020 £275		6" Aluminium front wheel 36 mm Froglegs Softroll (Only in silver) (Only with SH-f 52 cm)	
966401000??? £315		4" Aluminum front wheel 35 mm width Froglegs Softwheel (Not in stock - Longer delivery time)	Gold ☐ Purple ☐ Orange ☐ Green ☐
966401250??? £315		5" Aluminum front wheel 35 mm width Froglegs Softwheel (Not in stock - Longer delivery time)	


*** NOTE**

SH-f and SH-r are guidelines of heights, which we aim to meet. Though, these heights of the finished wheelchair may differ, depending on camber angle and choice of tires.

! Seat height at the rear SH-r (mandatory)

		Seat height at the rear (SH-r)								
Rear Wheel		39	40,5	41	43	44	45,5	46,5	48	
22"							-	-	-	 SH-r
24"		-	-					-	-	
25"		-	-	-						

Spoke Wheel							
	22"	24"	25"				
85102??01009	☐	☐	☐	£0	Spoke wheels (standard)		
8510??101002	☐	☐	☐	£215	Infinity Super Light 24 spokes (user weight limit 80 kg)		
85352??101008	-	☐	☐	£745	Black Spinergy 12 Spoke LX + Silver Hub		
85352??101208		☐	☐	£745	Blue Spinergy 12 Spoke LX+ Black Hub		
85352??101508		☐	☐	£745	White Spinergy 12 Spoke LX+ Black Hub		
85352??101308		☐	☐	£745	Red Spinergy 12 Spoke LX+ Black Hub		
85352??102108		☐	☐	£745	Black Spinergy 12 Spoke LX + Black Hub		
85352??101009	-	☐	☐	£580	Black Spinergy 18 Spoke SPOX+ Silver Hub		
85352??101209		☐	☐	£580	Blue Spinergy 18 Spoke SPOX + Black Hub		
85352??101509		☐	☐	£580	White Spinergy 18 Spoke SPOX + Black Hub		
85352??101309		☐	☐	£580	Red Spinergy 18 Spoke SPOX + Black Hub		
85352??101709		☐	☐	£580	Black Spinergy 18 Spoke SPOX+ Black Hub		
Wheels to order, extended delivery time, OBS color availability may change from manufacture. Please contact Wolturnus for more options.							
85352??101608	-	☐	☐	£795	Orange Spinergy 12 Spoke LX + Black Hub		
85352??101408		☐	☐	£795	Yellow Spinergy 12 Spoke LX + Black Hub		
85352??101808		☐	☐	£795	Pink Spinergy 12 Spoke LX + Black Hub		
85352??101708		☐	☐	£795	Green Spinergy 12 Spoke LX + Black Hub		
Tires							
	22"	24"	25"				
9662????00??	☐	☐	☐	£0	Std. Schwalbe tires, grey, 1"		
96622?01000?	☐	☐	☐	£95	Schwalbe Marathon Plus EVO, black 1"		
96622?010009	☐	☐	☐	£125	Solid tire, grey 1"		
96622?010021	-	☐	☐	£125	Solid tire, black 1"		
Pushrims							
85512?10200?	☐	☐	☐	£0	Aluminium pushrims Nature (standard)	Pushrims mounted short (std.)	☐
85512?102008	-	☐	☐	£160	Aluminium pushrims Black Ø19 6L		
85502?103009	-	☐	☐	£150	Stainless steel pushrims		
85552?000009	☐	☐	☐	£420	Titanium drivings		
8554202?0601	-	☐	☐	£265	Friction pushrims		
8554262?0601	-	☐	☐	£470	Curve L 6L Carbolife		
8554272?0601	-	☐	☐	£470	Curve L TG 6L Carbolife	Pushrims mounted long	☐
8554282?0601	☐	☐	☐	£570	Gekko 6L Carbolife		
9669282?0602	☐	☐	☐	£570	Gekko S 6L Carbolife		
8554302?0601	☐	☐	☐	£570	Quadro 6L Carbolife		
8554502?060?	-	☐	☐	£580	Oval Wood S BIG Light		
8554502?060?	-	☐	☐	£580	Oval Wood S BIG Dark		
85712?102009	☐	☐	☐	£150	Spoke guards transparent		



! Camber rear axle (Mandatory)		
611500010001	☐	1° Camber
611500030001	☐	3° Camber



! Rear Axle (Mandatory)									
		Adjustable rear axle - (Std. 8 cm) Pick RA below							
611000000001	!	RA:	☐ 2 cm	☐ 3 cm	☐ 4 cm	☐ 5 cm	☐ 6 cm	☐ 7 cm	☐ 8 cm

! Backrest System (Mandatory)	
!	BH: Back height delivered (BH 33-45 cm, adjustable) _____ cm



825010000000	☐	Straight back (std.)
82501000???? £465	☐	Wing Back HAWK ILSA: Integrated multi-adjustable back support Z: 3 cm
	☐	Wing Back HAWK ILSA: Integrated multi-adjustable back support Z: 5 cm
	☐	Wing Back HAWK ILSA: Integrated multi-adjustable back support Z: 7 cm


Anti-tipper (Recommended)		
953599992001 £275	☐	Swing away anti-tip - Right
953599992002 £275	☐	Swing away anti-tip - Left
		(Note that crutch holder can't fit together with anti-tip on small frames and 35/37,5 cm SW)




! Upholstery (mandatory)				
750000009000 £220	☐	W-fit seat cushion - Taylor made seat cushion for a perfect fit		!
- £0	☐	3D Knitted cover (Black, Std.)	Width / Depth	_____ cm / _____ cm
750000000005 £60	☐	Incontinens cover	Thickness	☐ 5cm Std. ☐ 7cm
			Firmness	☐ Medium Std. ☐ Hard


! Angle Adjustable Footrest (Mandatory)							
614500??0001	☐	VSL: Front of seat upholstery to rear edge of footrest Small: 35-46 cm Medium 36-47 cm Large 37-48 cm _____ cm				 USL/ VSL FRH-1 FRH-2	
	☐	FRH: Footrest height (Alternative to VSL)		FRH-1 _____ cm FRH-2 _____ cm			
95359990?70 £0	☐	Angle adjustable Hawk footrest (Note: FRD 14 cm deep) (std.)				 FRD: _____ cm	
829999010003 £435	☐	Footplate that can be flipped up to the rear (Fill out FRD)					
829999010009 £410	☐	Open footrest (Note: FRD: 8-15 cm) (Expect longer delivery time)				FRD: _____ cm	
840710001010 £705	☐	Freewheel 					
		FRA: Tilt on footrest in cm from the front to the back				 FRA	
		0 cm ☐	1 cm ☐	1,5 cm ☐	2 cm ☐		2,5 cm ☐

Sideguards (not optional with armrest)

Attachment system		Lift-off	Foldable on seat
			
Sideguard		95359930005C	681000010010
88100010000?	Dibond (std.)	☐ £0	☐ £115
88100020000?	Carbon	☐ £600	☐ £705
88100030000?	Carbon w. mudguard	☐ £655	☐ £765
767000000010	Thermal protection (3 mm)	☐ £145	

Armrest with ABS sideguards				
951000080220 88110010000? £460	☐	ABS-sideguards - 23,5 cm Armrest pad, short	 AS	
951000080126 88110010000? £460	☐	ABS-sideguards - 23,5 cm Armrest pad, staggered 2 cm		
951000080121 88110010000? £460	☐	ABS-sideguards - 27,5 cm Armrest pad		
951000080122 88110010000? £460	☐	ABS-sideguards - 27 cm Armrest pad, staggered 2 cm		
	!	AS: Height from seat tube to top of armrest (Note: e.g., With 24" wheels and SH-r 41-43 will min. be 23 cm and max. be 29 cm)		AS: _____ cm

Pushhandles			☐ NEIN / NEJ
885200011007 £435	☐	Height adjustable push handles	
885200051007 £480	☐	Foldable push handles mounted on consol with offset from back tubes	
8852000510?? £155	☐	Foldable push handles mounted back tubes (Does not work with ILSA back)	

Accessories			
953500003000 £150	☐	Double crutch holder (2 x Ø47 mm) Right	
	☐	Double crutch holder (2 x Ø47 mm) Left	
		(Note that crutch holder can't fit together with anti-tip on small frames and 35/37,5 cm SW)	
7600000100?? £60 £55	☐	Frame protectors front frame	
	☐	Frame protector crosstube on backrest	
76100000000? £135 £190	☐	Bag for catheter etc. (can be fixed under the seat) Small	
	☐	Bag for catheter etc. (can be fixed under the seat) Large	
7711000007?? £95 £95 £95	☐	BodyPoint calf strap Small 32" (81cm)	
	☐	BodyPoint calf strap Medium 36" (91cm)	
	☐	BodyPoint calf strap Large 40" (101cm)	
771100000747 £95	☐	BodyPoint Mobility Bag (medium)	

Bremsers (anbefalet valg)			☐ NEIN / NEJ
885110101001 £0	☒	Push-action wheel locks (Aluminium)	 
885110101003 £0	☐	Push-action wheel locks (Aluminium), bended handles (horizontal)	 
885110122002 £50 885110122002 £100	☐ ☐ ☐	Right extension handles for brake (Only for push-action wheel locks) Left extension handles for brake (Only for push-action wheel locks) Both sides - extension handles for brakes	
885110201009 £130	☐	Scissor wheel-locks with spring	 
885110200909 £165	☐	Scissor wheel-locks without spring	
885110101101 £165	☐	Push-action wheel locks (fibre plastic)	
885110201209 £190	☐	Scissor wheel locks small (fibre plastic)	

Subject to price changes, due to special combinations and checkbox errors	
Estimated Price (tax not included)	



WOLTURNUS

Kundeservice/Customer Service

Wolturnus A/S
Skalhuse 31
9240 Nibe

Phone/Telefon: +45 96 71 71 70

www.wolturnus.dk

Mail: ordre@wolturnus.dk

Produktion/Production

Wolturnus A/S
Halkjærvej 24B
9240 Nibe