



WOLTURNUS

Tukan

Rigid wheelchair
with open frame



»» I've always loved this minimalistic model. When you become one with your wheelchair, the possibilities are endless."

Birgitte Pedersen, wheelchair-user

Elegant and robust for an independent life

Tukan is a strong, rigid wheelchair with open frame that gives the wheelchair a light and elegant look. The open frame makes the wheelchair easy to clean and transport.

The Tukan wheelchair is a strong and durable rigid wheelchair, with a very robust frame as well as an elegant look. The frame is made of aluminium 7020, which is the strongest aluminium alloy that can be welded.

The open frame makes the Tukan easy to clean and to stow in a small car. That makes the Tukan a perfect choice for users, who live an active, mobile and independent life.

The Tukan series includes one standard and two customised models. The Tukan A Standard is designed to meet life's changing challenges. It's fully adjustable, so you can continually adapt the seating position and the wheelchair set-up.

The Tukan and the Tukan A are customised models, made to fit your measurements and meet your requirements. The Tukan has a fully welded backrest and rear axle, whereas the Tukan A comes with optional adjustability. This ensures a maximum flexibility regarding the individual adjustment of the seating posture.

Tukan key features

- › Active rigid wheelchair with open frame
- › Standard and customised models
- › Adapts to your requirements
- › Easy to handle and manoeuvre with minimum effort
- › Easy to clean and transport
- › Extremely high stability and low weight
- › Frames made of 7020 aluminium
- › Danish design and manufacture
- › Propulsion weight from 9 kg
- › User weight max 120 kg



› TUKAN

Customised for experienced users. Robust with fixed back and rear axle. Propulsion weight from 9 kg.



› TUKAN A STANDARD

Standard sizes with a range of seat widths. Fully adjustable. Propulsion weight from 9,5 kg.



› TUKAN A

Customised model with adjustable rear axle. Propulsion weight from 10 kg.