



WOLTURNUS

Tukan

Rigid wheelchair
with open frame



»» I've always loved this minimalistic model. When you become one with your wheelchair, the possibilities are endless."

Birgitte Pedersen, wheelchair-user

Elegant and robust for an independent life

Tukan is a strong, rigid wheelchair with open frame that gives the wheelchair a light and elegant look. The open frame makes the wheelchair easy to clean and transport.

The Tukan wheelchair is a strong and durable rigid wheelchair, with a very robust frame as well as an elegant look. The frame is made of aluminium 7020, which is the strongest aluminium alloy that can be welded.

The open frame makes the Tukan easy to clean and to stow in a small car. That makes the Tukan a perfect choice for users, who live an active, mobile and independent life. The Tukan and the Tukan A are customised models, made to fit your measurements and meet your requirements.

The Tukan has a fully welded backrest and rear axle. This makes it preferably suitable for experienced customers, who know all about their specific requirements and needs.

The Tukan A comes with optional adjustability. This ensures a maximum flexibility regarding the individual adjustment of the seating posture. The rear axle as well as the folding back and the footrest can be adjusted continuously.

Features

- › Active rigid wheelchair with open frame
- › Available with fixed back or folding back
- › Adapts to your requirements
- › Easy to handle and manoeuvre with minimum effort
- › Easy to clean and transport
- › Extremely high stability and low weight
- › Frame made of 7020 aluminium
- › Danish design and manufacture
- › Propulsion weight from 9 kg
- › User weight max 120 kg



› TUKAN

Customized for experienced users. Robust with fixed back and rear axle. Propulsion weight from 9 kg.



› TUKAN A

Tukan model with optional adjustability: adjustable rear axle, footrest & folding back adjustable by angle and height.



› OPEN FRAME

The open frame allows an easy storage in a car.