

W5, Tukan, Dalton Order Form

User information	
User Alias	
Supplementary Information	
Weight	Kg.
Height	cm.
Offer Only	☐
Order	☐
Payment Details	
EAN number	
Your PO number	
Invoice recipient	
Phone/Telefon	
E-Mail	
Contact Person	
Name and address	
Phone	
E-Mail	
Delivery Information	
Consultant	☐
User	☐
Dealer	☐
Other address	☐
Consultant / Date	

Project No.
Frame No.

Filled out by
Wolturnus

Custom-built wheelchair for active use Default Equipment

- Tailor-made rigid-frame wheelchair
- Fully welded or adjustable versions available
- Made of 7020 T6 aluminium
- Adjustable seat and back upholstery
- Impact-resistant anodized frame
- Rear wheels with quick-release axles


Comments:



Orderform is only compatible with
Adobe Acrobat Reader which is available at:

PC & Macbook
<https://get.adobe.com/uk/reader/>

Tablet & iPad
<https://www.adobe.com/acrobat/mobile.html>



Tailor-made Wheelchair

- ☐ This combination/option is available
- This combination/option is **not available**
- ! This information is mandatory and must always be filled out
-  Image Web-Link of Solution

For options or solutions not mentioned in this order form, please, contact Wolturnus directly.
See contact information at the bottom of this page.




* OBS! Tukan is not approved for front-frame mounted-handbikes and/or other electric devices mounted on the front frame of the wheelchair


Standard Frames							Max Userweight
11100000100?	£3,995	☒	1.	W5	Original (std. Ø25 seat tubes)		150kg max 48cm SW
11100000300?	£3,995	☐	2.	W5 S	Streamlined lower frame (std. Ø25 seat tubes)		150kg max 48cm SW
11100000400?	£3,995	☐	3.	W5 K	Short frame, the footrest is between casters (std. Ø25 seat tubes)		150kg max 48cm SW
11100000200?	£3,995	☐	4.	W5 D	Divided footrests (std. Ø30 seat tubes)		150kg max 48cm SW
11100000550?	£4,695	☐	5.	W5 XXL	Extra strong (std. Ø30 seat tubes) (For users above 150kg, SW > 48 cm)		250kg
120000001500	£3,995	☐	6.	Tukan A*	Adjustable rear axle frame (std. Ø30 seat tubes)		120kg max 45cm SW
120000002000	£3,995	☐	7.	Tukan B*	Welded rear axle shape B (std. Ø30 seat tubes)		120kg max 45cm SW
120000003000	£3,995	☐	8.	Tukan C*	Welded rear axle shape C (std. Ø30 seat tubes)		120kg max 45cm SW
120000004000	£3,995	☐	9.	Tukan D*	Welded rear axle shape D (std. Ø30 seat tubes)		120kg max 45cm SW
* OBS! Not approved for front-frame mounted- handbikes and/or other electric devices mounted on the front frame of the wheelchair							
14100000150?	£3,995	☐	10.	Dalton	Swing-away front frame (std. Ø25 seat tubes / front frame Ø30)		150kg max 48cm SW
141000001800	£4,695	☐	-	Dalton XXL	Swing-away front frame (std. Ø25 seat tubes / front frame Ø30)		Std. 160kg HD 250kg
+829999030030	£1,315	☐	11.	HD	Heavy Duty Swing-away front frame (std. Ø30 seat tubes footrest incl.)		250kg
Comments							

! Color (Mandatory)			
619000001009	£0	☐	Black anodized <i>no. 10</i>
619000001109	£0	☐	Silver anodized <i>no. 11</i>
619000001209	£0	☐	Blue anodized <i>no. 12</i>
619000001309	£0	☐	Red anodized <i>no. 13</i>
619000001409	£0	☐	Green anodized <i>no. 14</i>
619000001509	£0	☐	Bronze anodized <i>no. 15</i>
619000001609	£0	☐	Dove blue powder coat <i>no. 16</i>
619000001709	£0	☐	Antique green powder coat <i>no. 17</i>
619000001809	£0	☐	Blackberry powder coat <i>no. 18</i>
619000001909	£0	☐	Violet powder coat <i>no. 19</i>
619000002009	£0	☐	Steel blue powder coat <i>no. 20</i>
619000002109	£0	☐	White powder coat <i>no. 21</i>
619000002509	£0	☐	Shiny black powder coat <i>no. 25</i>
619000002609	£0	☐	Orange powder coat <i>no. 26</i>
619000002809	£0	☐	Matt black powder coat <i>no. 28</i>
619000003309	£0	☐	Dark Grey powder coat <i>no. 33</i>
619000003409	£0	☐	Ruby Red powder coat <i>no. 34</i>
619000003509	£0	☐	Brilliant Blue powder coat <i>no. 35</i>
619000003609	£0	☐	Aged Copper powder coat <i>no. 36</i>
619000003009	£625	☐	Custom colour (Powder coating) ! Colour code
	£0	☐	DECAL(S): The user's name or other message(s) printed on the frame
000000001102	£0	☐	1st. drawing of frame

OBS! Only the frame will be colored with selected paint/treatment. Other loose components, such as folding backrests, sideguards, armrest, suspension axle, etc. will be anodised black



Frame Data				
		!	SW: Seat width - outer seat tube measurement	SW: _____ cm
		!	FD: Frame depth	FD: _____ cm std. 44cm
		!	SD: Requested seat depth	SD: _____ cm std. 41cm
		☐	SD_{min}: Minimum seat depth (std. optimized in production)	SD_{min}: _____ cm optimized
		!	SH-r: Seat height from the rear edge of the seat to the floor	SH-r: _____ cm std. 42cm
		!	SH-f: Seat height from front of delivered seat depth to the floor	SH-f: _____ cm std. 48cm
		!	BH: Delivered Backrest Height (std. 35-41cm)	BH: _____ cm std. 35cm
		☒	Height adjustable backrest (std. +/-3cm, with 1,5cm increments)	+3cm / -3cm
		☐	Fixed backrest height	
 69100001010?		!	C: Camber angle rear axle 1°-3°-5°	C: _____ ° std. 3°
	£605	☐	RA: Fixed rear axle position	RA: _____ cm std. 8cm
		☐	RA: Adjustable rear axle position (Std. Frame = RA 3-13cm) (W5S-Frame = RA 6-9cm)	
		☐	HRA: Extra rear axle position (for handbike)	HRA: _____ cm std. 5cm
610000001009	£0	☐	Straight seat tubes	
610000??0001	£0	☐	ES: Ergonomic seat	ES: _____ cm std. 12cm
			BF: From the front of backrest to the front of footrest tube	BF: _____ cm std. 60cm
			BC: From backrest to the center of the caster house Std. BF-10cm Dalton FD+8cm	BC: _____ cm std. 50cm
			CW: Distance between the casters houses (FF+15cm, 4"castors)	CW: _____ cm FF+15cm
			TW: Maximum total width including wheels (SW+17cm, C=1°)	TW: _____ cm SW+17cm
612000??0009		☐	FB: U-V front frame (std.)	FB: _____ cm std. 15cm
613000??0009	!	☐	DF: U front frame, Length of seat tube up to bend (Max. FD - 9 cm)	DF: _____ cm std. 30cm
613000000009		☐	U front frame, straight down tube	
			FF: Front frame width Std. Ø25: 25,5 - 27,5 - 29,5 - 32,5 - 35,5 cm Std. Ø30: 25,0 - 27,0 - 29,0 - 32,0 - 35,0 cm	FF: _____ cm
Comments				




Backrest Data					
		!	BA: Backrest Angle (+/-)		BA: std. 0cm _____ cm
813100050???	£0	☒	1	Fully Welded Backrest	
813100051???	£305	☐	2	Folding Backrest	
691000010130	£60	☐	2+	Addon - Locking bracket for folding backrest	
813100053???	£415	☐	3	Light Style Folding Backrest (Max userweight 80kg/ Max SW 40cm)	
	£0	☒	Standard Backrest (No ILSA)		
8122002001??	£465	☐	Wing Back ILSA: Integrated Multi-adjustable (Std. Ø25)		Z = Depth Z: 0cm 3cm 5cm 7cm 10cm
8122002001??	£600	☐	Wing Back ILSA: Ball Clamped (fill out ZH, std. BH)(min. Z: 3cm)		ZH: std. BH _____ cm
			ZA: ILSA start bend		ZA: std. 19cm _____ cm
8122002001??	£770	☐	Custom ILSA Z: Custom ILSA Depth (Std. Ø25)		Z = Depth Z: std. 5cm _____ cm
	£770	☐	Custom Wing Back ILSA: Ball Clamp (Std. Ø25)		ZH: std. BH _____ cm
			ZA: ILSA start bend		ZA: std. 19cm _____ cm
825000000038	£0	☐	Default back brace shape		Depth: BD: std. 12cm _____ cm
825000000048		☐	Round back brace shape		Height: TH: std. 12cm _____ cm
812200200002	£205	☐	Narrow back (The Narrow back is not available in combination Integrated ILSA)		UBW: SW +/- 3cm _____ cm
812200200001		☐	Wide back (The wide back is not available in combination with the folding back or Integrated ILSA)		LBH: std. 20cm _____ cm
812200200020	£125	☐	a1: Bend on backrest Backwards (+) / Forwards (-)		a1: std. +1cm _____ cm
			a2: Bend on backrest Outwards (+) / Inwards (-)		a2: std. +0cm _____ cm
			b: Length of the upper, bent part of the back tubes		b: std. 15cm _____ cm
			c: Distance from the seat tube to bend		c: std. 20cm _____ cm



Comments




Footrest					
		☐	VSL: Front of seat upholstery to rear edge of footrest (std. 39-45 cm)	VSL: std. 42cm	cm
		☐	FRH: Height of the footrest (Alternative to VSL/LLL)	FRA: std. 1cm	cm
				FRH-1: std. 10cm	cm
				FRH-2: std. 9cm	cm
		☐	Height adjustable footrest (1,5cm increments Down+ / Up-)	+/-1,5	☐
		☐	Footrest height, Fixed by welding	+/-3	☐
8200FF001200	£0	☐	No Footrest or Custom Request		
8200FF001600	£0	☐	Angle adjustable footplate FRD 12cm (std.)	FRP:1cm	☐
82??99000000	£120	☐	Angle adjustable footplate FRD 16cm	FRP:3cm	☐
829999012001	£260	☐	Angle adjustable footrest Custom FRD (max. FRD 20 cm)	FRP:5cm	☐
829999012001	£295	☐	Angle adjustable footrest Custom FRD with bend (max. FRD 20 cm)	FRD: std. 12cm	cm
8299990000?9	£0	☐	2. Footrest, open with two tubes	HE: std. 1cm	cm
8299990100?1	£230	☐	3. Footplate welded on footrest		
8299990100?5	£190	☐	4. Footplate screwed on footrest w/tabs	FRD: std. 12cm	cm
8299990100?1	£380	☐	4. Carbon footplate screwed on footrest w/tabs		
8299990100?3	£435	☐	5. Footplate that can be flipped up to the rear (FRD min 15 cm)		
829999200005	£0	☐	6. Plastic Double footplates flip-up (Std. Dalton, Max 160kg)	FRD = 13cm	
9512000?000?	£280	☐	7. Alu Double footplates (Std. DaltonXXL HeavyDuty, max 250kg) (OBS! Max 180kg for flip-up version, above 180kg not flip-up)	FRD = 14cm	
835050001001	£810	☐	8. Elevating Footrest (Right Side) (Dalton model Exclusive Option)	FRD = 14cm	
	£810	☐	8. Elevating Footrest (Left Side) (Dalton model Exclusive Option)		
829999012001	£30	☐	HE: Height of bent edge on the footplate (only footplate 3, 4, 5)	HE: std. +1cm	cm
	£0	☐	FRP: Forward protrusion of footplate (only for option 2, 3, 4, 5 cm)	FRP: std. 0cm	cm
835099210001	60	☐	Rubber on footrest		
691000010001	60	☐	RH: Welded rings for leg strap, Height from floor.	RH: std. 23cm	cm
825000000003	£105	☐	Reinforced footrest for FreeWheel (only available footplate option 2, 3, 4)	FRD: 8-15 cm Max FRH 12cm & Max FRA 2cm	
8407100010??	£705	☐	Freewheel	FreeWheel Wolturnus Assembly	

Comments

Sideguards						
			Lift-off	Foldable OBS! (SW-s = SW+8mm)	Welded Sideguards OBS! Only on welded backrest	Screwed on Welded backrest OBS! In combination with Light Style Backrest SW-s = SW+14mm OBS! Not optional for Folding backrest
881000600003	£0	Standard ABS Sideguards (Std. with armrest)(Fixed shape)	☐			
88100010000?	£0	Custom Dibond Sideguards (Std.)	☐	☐		☐
881000?0000?	£280	Plastic Sideguards w/ Mudguards (MC not optional)	☐	☐		☐
88100040000?	£250	Aluminium Sideguards	☐	☐	☐	☐
8810002000??	£600	Carbon Sideguards	☐	☐		☐
88100050000?	£620	Aluminium Sideguards with Mudguards	☐	☐	☐	
8810005000022	£620	With Black sideguards adapter screw				☐
8810005000012	£620	With Nature sideguards adapter screw				☐
8810005000032	£620	With Red sideguards adapter screw				☐
8810005000042	£620	With Blue sideguards adapter screw				☐
8810005000052	£620	With Green sideguards adapter screw				☐
8810005000062	£620	With Bronze sideguards adapter screw				☐
8810005000005	£620	With Blank Light Style Backrest screws				☐
883124101034	£870	In Combination with armrest	☐		☐	☐
8810003000??	£655	Carbon Sideguards with Mudguards	☐	☐		
8810003000022	£655	With Black sideguards adapter screw				☐
8810003000012	£655	With Nature sideguards adapter screw				☐
8810003000032	£655	With Red sideguards adapter screw				☐
8810003000042	£655	With Blue sideguards adapter screw				☐
8810003000052	£655	With Green sideguards adapter screw				☐
8810003000062	£655	With Bronze sideguards adapter screw				☐
8810003000005	£655	With Blank LS sideguards adapter screws				☐
883124101030	£845	In Combination with armrest	☐			☐
-	£0	No Sideguards			☐	
681000010007	£0	☐ Addons for Lift-off Solution:				
681000010004	£65	☐ Bracket Fixed with Screw				
681000010004	£150	☐ 1,5cm Expansion kit, adds 3cm to SW				
883124101014	£200	☐ Double pin, extra balusterholder (For reinforcing removable sideguards/armrest).				
883124101014	£200	☐ Adjustable console for lift-off sideguards (Not in combination with armrest)				
681000010052	£65	☐ Addons for welded-/foldable-/screwed on sideguards:				
681000010052	£65	☐ Reinforced with extra tab on seat				
		GH: Height from tire to underside of sideguards (clearance at 8 bar) GL: Length from backrest to front of sideguards (std. follow tire) MC: Mudguard cover (std. center of tire)(max 4cm sideguards)	GH: std. 1,0cm _____ cm GL: follow tire _____ cm MC: center _____ cm			
767000000010	£145	☐ 3mm Thermal protection (on both right and left sideguard)	SW-s = SW-6mm			
767000000013	£145	☐ 10mm Thermal protection (on both right and left sideguard)	SW-s = SW-20mm			
767000000014	£145	☐ 20mm Thermal protection (on both right and left sideguard)	SW-s = SW-40mm			
Comments						






Armrest (incl. std. ABS sideguards)					
	£0	☐	-	No Armrest	
881001000003	£460	☐	1.	23,5 cm Armrest pad, short	AB = 29 cm (Max. 150kg)
881001100003	£460	☐	2.	27,5 cm Armrest pad	AB = 30,5 cm (Max. 150kg)
881001200003	£460	☐	3.	23,5 cm Armrest pad, displaced by 2 cm	AB = 27,5 cm (Max. 150kg)
881001300003	£460	☐	4.	27 cm Armrest pad, displaced by 2 cm	AB = 30 cm (Max. 150kg)
881001700003	£665	☐	5.	37 cm Armrest pad, long, displaced 2 cm	AB = 41 cm (Max. 150kg)
881001700003	£790	☐	6.	Reinforced 37cm armrest pad, long, displaced 2 cm	AB = 41 cm (Max. 250kg)
881001700003	£705	☐	7.	26 cm leather armrest pad, width 8 cm	AB = 31 cm (Max. 150kg)
881001700003	£830	☐	8.	Reinforced 26 cm leather armrest pad, width 8 cm	AB = 31 cm (Max. 250kg)
881001700003	£705	☐	9.	31 cm armrest pad, leather, width 8 cm	AB = 34,5 cm (Max. 150kg)
881001700003	£830	☐	10.	Reinforced 31 cm armrest pad, leather, width 8 cm	AB = 34,5 cm (Max. 250kg)
883124101002	£465	☐	11.	Swing-away armrest, attached behind back (fill out AB)	AB: _____ cm std. 35cm
AS: Armrest to seat (armrest height)					AS: _____ cm std. 26cm




Pushhandles					
81310005??0?	£0	☐	-	No Push handles	
81310005?106	£260	☐	1.	Push handles fixed by screws (fixed backrest height)	
81310005?206	£310	☐	1.	Push handles fixed by screws (adjustable backrest height)	
81310005?108	£260	☐	2.	Push handles welded on (fixed backrest height)	
81310005?208	£310	☐	2.	Push handles (adjustable backrest height)	
81310005?107	£155	☐	3.	Foldable push handles, mounted in top of back (fixed backrest height)	
81310005?207	£195	☐	3.	Foldable push handles, mounted in top of back (adjustable backrest height)	
81310005??02	£460	☐	4.	Adjustable push handles Dia. Ø22 mm	
81310005??03	£460	☐	4.	Adjustable push handles Dia. Ø25 mm (for lifting)	
81310005??04	£310	☐	5.	Push handles welded on extra tubes welded behind backrest	
81310005??05	£480	☐	6.	Foldable push handles, behind backrest	
PH: Height from floor to anticipated position (Only available for option 4, 5 et 6)					PH: _____ cm std. SH + BH


Comments


Wheel-locks					
	£0	☐	0.	No wheel-locks	         
885110101001	£0	☐	1.	Push activated wheel-locks (aluminium)	
885110101003	£0	☐	2.	Push-action wheel-locks (aluminium), bent handles (horizontal)	
885110111001	£165	☐	3.	Pull-activated wheel-locks (aluminium)	
885110122002	£50	☐	+	Extension handles for wheel-locks - Mounted on the Left - for Alu	
885110122002	£50	☐	+	Extension handles for wheel-locks - Mounted on the Right - for Alu	
885110201009	£130	☐	4.	Scissor wheel-locks with spring	
885110200909	£165	☐	5.	Scissor wheel-locks without spring	
885110101101	£165	☐	6.	Push-action wheel-locks (plasticfibre)	
885110111101	£165	☐	7.	Pull-activated wheel-locks (plasticfibre)	
951519204000	£75	☐	+	Extension handles for - Mounted on the Left - for Fiber	    
951519204000	£75	☐	+	Extension handles - Fiber- Mounted on the Right - for Fiber	
885110201209	£190	☐	8.	Sport Scissor Fiberplastic	
885110201300	£210	☐	9.	Wing Scissor Aluminium	
840000207242	£375	☐	10.	One-hand activated wheel-lock from OB (Right)	
840000207239		☐		One-hand activated wheel-lock from OB (Left)	
9515????00??	£0	☐	x	Mounted with clamps <i>(if no wheellock is selected, a loose set clamps will be sent)</i>	
691000020001	£145	☐	y	Extra tubes welded on frame to for wheel-locks <i>(Not optional for Dalton) / (included in WSSL)</i>	
-	£145	☐	z	Screwed on welded bushing <i>(only optional for Wing scissor)</i>	
	£0	☐	-	No Wheellock Mounting	

Anti-tip (Recommended) & Rear Axle Attachments					
	£0	☐		No Anti-tip	 Height from floor Distance to axle AF: 5-6cm _____ cm AW: 35-36cm _____ cm
8861101013??	£275	☐	a.	Anti-tip w/ Quick Release - Right <i>(Max 120kg)</i>	
	£275	☐		Anti-tip w/ Quick Release - Left <i>(Max 120kg)</i>	
	£550	☐		Anti-tip w/ QR - Right & Left <i>(Max 120kg)</i>	
	£275	☐		Anti-tip w/ QR - Center mounted <i>(Only for welded rear axle) (Max 120kg)</i>	
88611010130?	£275	☐	b.	Swing away anti-tip - Right <i>(Max 120kg)</i>	
	£275	☐		Swing away anti-tip - Left <i>(Max 120kg)</i>	
	£550	☐		Swing away anti-tip - Right & Left <i>(Max 120kg)</i>	
	£275	☐		Swing away anti-tip Center mounted <i>(Only for welded rear axle) (Max 120kg)</i>	
8861101013??	£355	☐	c.	Swing away anti-tip - Right Reinforced	
	£355	☐		Swing away anti-tip - Left Reinforced	
	£710	☐		Swing away anti-tip - Right & Left Reinforced	
691000020020	£230	☐		Tip assistant pedal Right Side <i>(Max 120kg)</i>	 
691000020020	£230	☐		Tip assistant pedal Left Side <i>(Max 120kg)</i>	
886210101001	£285	☐		Double crutch holder (2 x Ø47 mm) Right	
886210101001	£285	☐		Double crutch holder (2 x Ø47 mm) Left	

Comments




Neck Support						
	£0	☐	No head support		NH: _____ cm std. 32cm	
884100013001	£200	☐	Prepared for head support		HH: _____ cm std. 32cm	
88410001400?	£290	☐	Headsupport, adjustable - choose cushion below <i>(mandatory for crash test)</i>		HD: _____ cm std. 32cm	
951000080130	£80	☐	a.	Small rounded (23 x 10 x 4 cm)		
951000080131	£95	☐	a.	Large rounded (30 x 14 x 4 cm)		
951000080141	£95	☐	b.	Small square (14 x 8 x 5 cm)		
951000080142	£95	☐		Large square (22 x 10 x 5 cm)		
840000100831	£180	☐	c.	OB Comfort Plus, Str. S (20 x 14 x 6 cm)		
840000100832	£180	☐		OB Comfort Plus, Str. M (25 x 15 x 6 cm)		
884100250100	£160	☐	SF: Attachment for body strap fixation		SF: _____ cm std. 32cm	
Car Fixation & Reinforcement options						
825000000001	£240	☐	Extra reinforced backrest <i>(mandatory for crash test)(Not optional for Light Style Backrest)</i>			
825000000009	£335	☐	Folding Backrest Heavy Duty Reinforcement <i>(dual fold)</i> <i>(for XXL chairs with foldable backrest and crash test, or 150kg very active user or 170kg user)</i>			
88700020000?	£415	☐	4-point fittings for seat belt fastening <i>(Not optional for Tukan or Light Style Backrest)</i>			
887000200011	£1,995	☐	Fittings for Dahl Docking Station <i>(Not optional for Tukan or Light Style Backrest)</i>		Standard DH: _____ cm 6cm - 8,5cm - 10cm DH: _____ cm std. 8,5cm DL: _____ cm 23cm	
Comments						

Upholstery					
	£0	☐	1.	No Backrest Upholstery	
710000001110	£0	☐	2.	Backrest Comfort Upholstery (Std. 3D Mech Fabric with logo)	
-	£0		+	Extra wide back upholstery (backrest cover)	☐ +0 ☐ +1 ☐ +2 ☐ +3
-	£0		+	Extra wide backrest straps (for a deeper shape of back)	☐ +0 ☐ +1 ☐ +2 ☐ +3
710000001120	£260	☐	3.	Custom Backrest Upholstery (Std. 3D Mech Fabric with logo)	
710099990010	£365	☐		Custom Leather Backrest upholstery, (without logo)	
710000002502	£0		+	Extra wide back upholstery (backrest cover) (Note: +4 & +5cm, will have cutout fl crosstube)	☐ +0 ☐ +1 ☐ +2 ☐ +3 ☐ +4 ☐ +5
-	£0		+	Extra wide backrest straps (for a deeper shape of back)	☐ +0 ☐ +1 ☐ +2 ☐ +3
761000000010	£10	☐	+	Pocket on lower back upholstery (w/zipper)	
750000000103	+ £	☐	+	Comments (per request):	
710099990001	£0	☐	☒	Std. Backrest top w/ leather corners	
710000001?20	£165	☐	☒	Custom backrest top Comments (per request):	
710000000000	£0	☐	1.	Std. Seat Upholstery	
725029351009	£130	☐	2.	Dibond plate, on oval tubes (fixed with velcro, +/-10mm adjustable.)	
725029351009	£130	☐	3.	Dibond Seat fixed with Rivets (No Adjustability)	
-	£0	☐	4.	No Seat Upholstery or Dibond Seat	
Cushion					
750000009000	£220	☐	1.	W-fit seat cushion - Taylor made to fit seat	
750000000001	£0	☐	+	3D Knitted cover (Black, Std.)	
750000000005	£60	☐	+	Incontinens cover (Fabric side out)	Width: W: std. SW _____ cm
750000000006	£60	☐	+	Incontinens cover (Smooth side out)	Depth: D: std. SD _____ cm
		☒		Softness Medium Cold Foam (Select Thickness)	5cm ☐ 7cm ☐
		☐		Softness Hard Granular Foam (Select Thickness)	5cm ☐ 7,5cm ☐
750000000000?	£60	☐		Extra Cushion cover (same as ordered)	
Accessories					
760000010001	£60	☐	1.	Frame protectors front frame	
760000010015	£55	☐	2.	Frame protector crosstube on backrest	
761000000001	£135	☐	3.	Bag for catheter etc. (can be fixed under the seat) Small	
761000000002	£190	☐		Bag for catheter etc. (can be fixed under the seat) Large	
771100000751	£95	☐	4.	BodyPoint calf strap Small 32" (81cm)	
771100000750	£95	☐		BodyPoint calf strap Medium 36" (91cm)	
771100000749	£95	☐		BodyPoint calf strap Large 40" (101cm)	
771100000747	£95	☐		BodyPoint Mobility Bag (medium)	

		4"	5"	6"	Castors (*Extended delivery time)	
860010101401	£0			☐	Pneumatic Grey Ø140	
860010101253	£0		☐		Pneumatic Grey 40mm	
8600102009?0	£150	☐	☐		Aluminium	
8600101020???	£0	☐	☐		Black Fiber plastic 34 mm width	
	£0	☐			Blue Fiber plastic 34 mm width	
	£0	☐			Red Fiber plastic 34 mm width	
	£0	☐			White Fiber plastic 34 mm width	
	£0	☐			Silver Fiber plastic 34 mm width	
860020000000	£45		☐		Anti-haircap for fiber plast wheels	
860010201???	£275	☐	☐	☐	Silver Aluminum 36 mm Froglegs Softroll	
	£275	☐	☐		Black Aluminum 36 mm Froglegs Softroll	
	£275	☐	☐		Blue Aluminum 36 mm Froglegs Softroll	
	£275	☐	☐		Red Aluminum 36 mm Froglegs Softroll	
	£315	☐	☐		*Purple Aluminum 36 mm Froglegs Softroll	
	£315	☐	☐		*Orange Aluminum 36 mm Froglegs Softroll	
		+		Anti-haircap for froglegs are incl.		
860010202???	£275	☐	☐		Black Aluminum 36 mm Froglegs Whisper	
	£275	☐	☐		Silver Aluminum 36 mm Froglegs Whisper	
		+		Anti-haircap for froglegs are incl.		
Castor Forks						
	£0		☐		Black Aluminum Front Fork Color <i>(Std.)</i>	
	£135		☐		Silver Aluminum Front Fork	
	£135		☐		Red Aluminum Front Fork	
	£135		☐		Blue Aluminum Front Fork	
955067070012	£500		☐		Castors Pin Locks, (Cable activated)	
955067070011	£500		☐		Castors Pin Locks, (Lever activated)	
862400001009	£260		☐		Fork w. Quick Release Axle	
Comments						


		22"	24"	25"	26"	Tires	ERTRO
	£0	☐	☐	☐	☐	No Tires or Custom (Fill out ERTRO/Name)	
96622?0100??	£0	☐	☐	☐	☐	Schwalbe Right Run grå, 1" (Std.)	(25-489) (25-540) (23-559) (25-590)
96622?01000?	£95	☐	☐	☐	☐	Schwalbe Marathon Plus EVO, black 1"	(25-489) (25-540) (23-559) (25-590)
96622?010009	£125	☐	☐	☐	☐	Solid tire grey 1"	(25-489) (25-540) (23-559)
966224010021	£125		☐			Solid tire black 1"	(25-540)
966224010013	£60		☐			Schwalbe Right Run Plus, 1"	(25-540)
966224011008	£60		☐			Schwalbe Marathon black, 1-3/8"	(37-540)
966224010007	£60		☐			Schwalbe Marathon grey, 1-3/8"	(37-540)



		22"	24"	25"	26"	DriveWheels
	-£50	☐	☐	☐	☐	No Drivewheels or Custom (Fill out Name)
85102??01009	£0	☐	☐	☐	☐	Spoke wheels (standard)
85102??0100?	£0		☐	☐		Spoke wheel radial
8510??0101002	£215	☐	☐	☐		Infinity Super Light 24 spokes (user weight limit 80 kg)
851124201009	£395		☐			MTB Wheel 36 Spoke (incl. tire and pushrims)
8537281110??	£680	☐	☐		☐	Drum Brake w/ cable lock (Manual Lock)
	£680	☐	☐		☐	Drum Brake w/ cable lock (Automatic Lock)
8537281110??	£1,090	☐	☐			One Hand Drive (Right Side)
	£1,090	☐	☐			One Hand Drive (Left Side)

		22"	24"	25"	26"	Spinergy (*Extended delivery time)	
85352??01?09	£745		☐	☐	☐	Black Spinergy 12 Spoke LX + Silver Hub	
85352??012108	£745		☐	☐	☐	Black Spinergy 12 Spoke LX + Black Hub	
85352??0101208	£745		☐	☐	☐	Blue Spinergy 12 Spoke LX+ Black Hub	
85352??0101508	£745		☐	☐	☐	White Spinergy 12 Spoke LX+ Black Hub	
85352??0101308	£745		☐	☐	☐	Red Spinergy 12 Spoke LX+ Black Hub	
85352??0101608	£795		☐	☐	☐	*Orange Spinergy 12 Spoke LX + Black Hub	
85352??0101408	£795		☐	☐	☐	*Yellow Spinergy 12 Spoke LX + Black Hub	
85352??0101808	£795		☐	☐	☐	*Pink Spinergy 12 Spoke LX + Black Hub	
85352??0101708	£795		☐	☐	☐	*Green Spinergy 12 Spoke LX + Black Hub	
85352??0101009	£580		☐	☐	☐	Black Spinergy 18 Spoke SPOX+ Silver Hub	
85352??0101?09	£580		☐	☐	☐	Black Spinergy 18 Spoke SPOX+ Black Hub	
85352??0101209	£580		☐	☐	☐	Blue Spinergy 18 Spoke SPOX + Black Hub	
85352??0101509	£580		☐	☐	☐	White Spinergy 18 Spoke SPOX + Black Hub	
85352??0101309	£580		☐	☐	☐	Red Spinergy 18 Spoke SPOX + Black Hub	
85362??241008	£795		☐	☐	☐	*Black Spinergy 18 Spoke SLX18+ Silver Hub	
85362??181208	£795		☐	☐	☐	*Blue Spinergy 18 Spoke SLX18 + Black Hub	
85362??181508	£795		☐	☐	☐	*White Spinergy 18 Spoke SLX18 + Black Hub	
85362??181308	£795		☐	☐	☐	*Red Spinergy 18 Spoke SLX18 + Black Hub	
85352??181107	£795		☐	☐		*Spinergy Black Hub - B-LXL	

DriveWheel Accessory						
953001040120	£90		☐			Quick-Release axle with Tetra Grip

Comments

		22"	24"	25"	26"	Pushrims	
	£0	☐	☐	☐	☐	No Pushrims or Custom	
85512?10200?	£0	☐	☐	☐	☐	Aluminium Nature <i>(std.)</i>	
85502?103009	£150		☐	☐	☐	Stainless steel pushrims	
85552?000009	£420	☐	☐	☐	☐	Titanium drivring	
8554202?0601	£265	☐	☐	☐	☐	Friction pushrims	
8554262?0601	£470		☐	☐	☐	Curve L Carbolife 6L	
8554272?0601	£470		☐	☐	☐	Curve L - TG Carbolife 6L	
8554282?0601	£570	☐	☐	☐	☐	Gekko Carbolife Black	
9669282?0602	£570	☐	☐	☐	☐	Gekko S Carbolife Black	
8554302?0601	£570	☐	☐	☐	☐	Quadro Carbolife	
8554502?060?	£580		☐	☐		Oval Wood S BIG DARK	
8554502?060?	£580		☐	☐		Oval Wood S BIG LIGHT	
					☐	Pushrim Mounted Short <i>(std.)</i>	
					☐	Pushrim Mounted Long	
85712?102009	£150		☐			Transparant Spoke Guards	


Electric Support (Prices doesnt incl. electric equipment)						
		i	x: Only prepared for customers equipment			
	£0	☐	No Electric Wheels			
853727110001	£1,030	☐	24" Decon E-Move			R: ☐ L: ☐ JB: std. 35cm _____ cm JA: std. 4cm _____ cm
853727110002	£1,030	☐	24" Decon E-Drive (Joystick)			
853727110004	£1,030	☐	24" Alber E-Motion M25			
853727110005	£1,030	☐	24" Alber Twion (Joystick)			
853727110007	£1,030	☐	24" Alber E-Fix (Joystick)			
951099250460	£170	☐	Bracket f/ Alber/Decon Anti-tip (Right Side)			
951099250560	£170	☐	Bracket f/ Alber/Decon Anti-tip (Left Side)			
	£0	☐	No Rear Axle Motors			
853728120000	£195	☐	Smart Drive			
	£195	☐	Smoov			
	£0	☐	No Handbikes			
853728111039	POA	☐	Other/Custom			
853728111062	£485	☐	Batec			
853728111061	£485	☐	Klaxon Klick			
853728111064	£475	☐	24" Decon Infinity			
853728111063	£475	☐	24" Decon Liberty			
825000000007	£95	☐	Reinforced Crosstube f/ handbike (Ø25) CTP: Reinforced Crosstube Placement CTD: Reinforced Crosstube Depth			CTP: std. 4cm _____ cm CTD: std. 5cm _____ cm

Comments

Subject to price changes, due to special combinations and checkbox errors	
Estimated Price (tax not included)	

Custom Solutions and Equipment

000000001104 £90	☐	MDR: Risk mangement analysis of custom- equipment and/or solution on wheelchair/medical device. (Required for custom solutions and/or combinations of equipment, made by, mounted by or prepared for by Wolturnus A/S). Note: Not all standard/MDR validated solutions are mentioned in this orderform, please contact Wolturnus A/S for further information. (Price of the custom solution is not included)	EN ISO 14971
P.4			
P.5			
P.6			
P.7			
P.8			
P.9			
P.10			
P.11			
P.12			
P.13			
P.14			



How to use Adobe Acrobat on Mobile Devices (Phones and Tablets):

How to send an orderform from your mobile device:

1. Launch the "Adobe Acrobat" app on your device
2. Find your document on your list and press the "three dots" icon :
3. Press the share icon  (several options is available)

How to rename a orderform on your mobile device:

1. Launch the "Adobe Acrobat" app on your device
2. Go to the "files"  tab at the bottom menu
3. Find your document and press the "three dots" icon :
4. Press Rename, A

How to make a copi of an orderform:

1. Launch the "Adobe Acrobat" app on your device
2. Find your document on your list and press the "three dots" icon :
3. Press "save copi" 

Note: Please avoid using the "print to pdf" feature available on PC