



WOLTURNUS

W5

Strong active
rigid wheelchair



» I've tried everything in a W5. I've travelled all over the world and tested it to the limit. The W5 is reliable, has a perfect fit, and it doesn't break.

Jens Holst, wheelchair-user

Strong and stable for an active life

The W5 is a light and robust rigid wheelchair, designed for an active life. It has the optimum strength-to-weight ratio and makes optimum use of your energy.

The W5 is a strong, rigid wheelchair designed for an active lifestyle. The frame is made of aluminium 7020 alloy, which is the strongest aluminium alloy that can be welded.

The W5 series includes a total of eight standard and customised models. The standard wheelchairs are fully adjustable and come in standard sizes. The customised W5 models are made to fit each individual wheelchair user's measurements and requirements.

The W5 SL has a very light frame as well as slim frame tubes, which gives this wheelchair an elegant look. The W5 XXL, on the other hand, comes with a reinforced frame. It is made for active users, whose high weight and physical strength ask for special requirements.

Visit our homepage to see the full range of our W5 series.

W5 key features

- › Active rigid wheelchair with boxed frame
- › Strong and robust for an active life
- › Standard and customised models
- › Adapts to your requirements
- › Extremely high stability and low weight
- › Easy to handle and manoeuvre with minimum effort
- › The frames are made of 7020 aluminium
- › Danish design and manufacture
- › Propulsion weight from 6 kg
- › User weight max 80-250 kg



› W5

Robust wheelchair with rigid back and welded rear axle for experienced users. Max. user weight 170 kg.



› W5 SL

Elegant, light frame with narrow frame tubes. Max. user weight 80 kg.



› W5 XXL

Reinforced model, for larger and physically strong users. Max. user weight 250 kg.