

YOGA



Seat width	36 - 39 - 42 - 45 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Front from 43 a 53 cm Rear from 41 a 49 cm (adjustable)
Footplate distance	From 35 to 54 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	90° - 100°
Total length	90 cm approx. (conf. std)
Total width	Seat width + 17 cm (with 0° camber)
Encumbrance	28 cm approx. (with 0° camber)
Backrest angle	4 angles (adjustable)
Camber	0° - 2°
Weight without rear wheels	Starting from 9,6 kg approx.
Wheelchair weight	Starting from 12,2 kg approx.
Maximum load	125 kg