

International TiLite Z Series Order Form

Effective January 7, 2022 (New Items in Red)

MONO-TUBE RIGID FRAME

Page 1 of 8

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:	User Weight:	Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Please Email to the Permobil Manual Customer Service Team at: TiLite_Global_Support@Permobil.com

FRAME STYLE

AZFS1	<input type="checkbox"/> Aero Z Aluminum Frame (Full Adjustability. 265 Pound Weight Limit. 150 Pound Weight Limit for Rear Seat Width A and Seat Depth B of 10"x10" to 12"x12". For Rear Seat Width A and Seat Depth B of 15" or greater, use K0005. For Rear Seat Width A and Seat Depth B of less than 15", use E1235.).....	£1,912
Z2FS2	<input type="checkbox"/> ZRA Titanium Frame (Full Adjustability. 265 Pound Weight Limit. 150 Pound Weight Limit for Rear Seat Width A and Seat Depth B of 10"x10" to 12"x12". For Rear Seat Width A and Seat Depth B of 15" or greater, use K0005. For Rear Seat Width A and Seat Depth B of less than 15", use E1235.).....	£2,705
Z2FS1	<input type="checkbox"/> ZR Titanium Frame Ultimate TiFIT (Minimal Adjustability. 265 Pound Weight Limit. 150 Pound Weight Limit for Rear Seat Width A and Seat Depth B of 10"x10" to 12"x12". For Rear Seat Width A and Seat Depth B of 15" or greater, use K0005. For Rear Seat Width A and Seat Depth B of less than 15", use E1235.).....	£2,977
ZHDR1	<input type="checkbox"/> Heavy Duty Frame (Not available on Aero Z. Available on ZR and ZRA only. 300 Pound Weight Limit. No charge when 20" Rear Seat Width A and 20" Seat Depth B are selected. Includes Tension Adjustable by Straps Back Upholstery. Requires Folding Aluminum Adjustable Height Backrest, Billet Aluminum Front Wheels, Shadow or Spinergy SPOX Rear Wheels. LX Rear Wheels only available with 0 or 2 degree camber. Not available with SuperLite Edition, Slipstream Fork, Velcro Adjustable Upholstery, Swing Away Tubular Armrest or Swing Away Flip Back Tubular Armrest.).....	£324
ZRFR1	<input type="checkbox"/> Reinforced Frame (No charge when 20" Rear Seat Width A and 20" B Seat Depth are selected.).....	£253

TRANSIT OPTION

(Please refer to the ZRA Transit Tie-Down Option Brochure for important information regarding the Transit Tie-Down Option. Not Available with ZR.)

ZTTDO3	<input type="checkbox"/> Transit Tie-Down (Consists of four Tie-Down Brackets, two Pelvic Belt Mounting Brackets and Reinforced Frame.).....	£257
ZWTB	<input type="checkbox"/> Transport Bracket Option (Consists of four Tie-Down Brackets. Unoccupied use only.).....	£257

FRAME FINISHES

ZFF1	<input type="checkbox"/> Satin (Available with Titanium Frame only.).....	STD	
ZSCBT1	<input type="checkbox"/> Satin Finish Buffing Pads.....	£11	
ZFF3	<input type="checkbox"/> Tattooed Titanium™ and Aluminum (Tattoos on front frame tubes. Titanium frame finished in Satin. Aluminum frame is tattooed on blasted aluminum with a clear coat. Tattoos will be trimmed and adjusted to fit your frame.).....	£257	
ZFF45	<input type="checkbox"/> Need for Speed <div style="display: inline-block; width: 200px; border-bottom: 1px solid black; margin-top: 5px;"> </div>	ZFF44 <input type="checkbox"/> Flower Garden <div style="display: inline-block; width: 200px; border-bottom: 1px solid black; margin-top: 5px;"> </div>	
ZFF43	<input type="checkbox"/> Mosaic <div style="display: inline-block; width: 300px; border-bottom: 1px solid black; margin-top: 5px;"> </div>		

FRAME FINISHES

High Gloss Paint	<input type="checkbox"/> ZFF8 High Gloss Black	<input type="checkbox"/> ZFF10 TiLite Red	<input type="checkbox"/> ZFF29 Neon Pink
Metallic Paint	<input type="checkbox"/> ZFF14 Ocean Blue Metallic	<input type="checkbox"/> ZFF15 Tangerine Metallic	<input type="checkbox"/> ZFF30 Acid Green Metallic
Pearlescent Paint	<input type="checkbox"/> ZFF17 Midnight Blue Pearl	<input type="checkbox"/> ZFF18 White Pearl	
Matte Paint	<input type="checkbox"/> ZFF19 Black Lava	<input type="checkbox"/> ZFF27 Titanium Grey	
Candy Paint	<input type="checkbox"/> ZFF22 Ruby Red	<input type="checkbox"/> ZFF23 Azure Blue	<input type="checkbox"/> ZFF25 Electric Plum
Matte Textured Paint	<input type="checkbox"/> ZFF40 Hyper Blue	<input type="checkbox"/> ZFF41 Flame Orange	<input type="checkbox"/> ZFF42 OD Green..... STD

COLOR ANODIZE PACKAGE

Color Anodize Package includes: **Caster Mount, Caster Fork, and Caster Caps**. If no anodize color is selected, the default will be black.)

<input type="checkbox"/> ZANO1 Black.....			STD
<input type="checkbox"/> ZANO12 Gold	<input type="checkbox"/> ZANO15 Green	<input type="checkbox"/> ZANO13 Purple	<input type="checkbox"/> ZANO14 Orange
<input type="checkbox"/> ZANO10 Blue	<input type="checkbox"/> ZANO11 Red	£40	



Permobil Ltd. Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN
 Email: Customer_support_uk@permobil.com
 Phone: +44 (0)1484 722 888
 VAT No: 780130457
 Company Reg. No: 4284728

International TiLite Z Series Order Form

TiFit: TAILORED INDIVIDUALIZED FRAME

REAR SEAT WIDTH A Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

- ZSW4 10" (25.4cm) 11" (27.9cm) 12" (30.5cm) 13" (33cm) 14" (35.6cm) 15" (38.1cm)
 16" (40.6cm) 17" (43.2cm) 18" (45.7cm) 19" (48.3cm),.....**STD**
 ZSW2 20" (Greater than 20" size avail. through CMD),.....**£368**

FRONT SEAT WIDTH (SEAT TAPER) T Measure from inside of front frame tube to inside of opposite front frame tube. (Required with V-Front End v. Not Available with Standard Front End I.)

- ZFSW1 No Taper (2.5" (6.4cm) Narrower than Rear Seat Width A.),.....**STD**
 ZFSW3 1.5" (3.8cm) Taper (4" (10.2cm) Narrower than Rear Seat Width A. Includes Tension Adjustable by Straps Seat Upholstery. Not available with 10" (25.54cm) and 11" (27.9cm) Seat Widths A.),.....**£209**

SEAT DEPTH B Measure from front of back post to the front edge of seat sling. Frame bend will start .5" in front of seat sling unless Custom Frame Depth is specified, or the combination of specifications affect the rolling dynamics of the chair.

- ZSD1 10" (25.4cm) 11" (27.9cm) 12" (30.5cm) 13" (33cm) 14" (35.6cm)
 15" (38.1cm) 16" (40.6cm) 17" (43.2cm) 18" (45.7cm) 19" (48.3cm),.....**STD**
 ZSD2 20" (50.8cm),.....**£388**
 ZSD3 Custom Frame Depth: _____" Frame Depth Longer or Shorter Than Seat Depth D Measure from front of seat sling to beginning of frame bend. Ranges from 0.5" to 2" and can be given as specific as 1/4" increments. A Custom Frame Depth outside specified range may be available through our Creative Mobility Division. A Custom Frame Depth of 20" will have an additional charge.),.....**NCO**

ERGONOMIC SEAT Measure from front of back post to beginning of upward bend of seat tube.

- ZERG1 ZR Ergonomic Seat (Not available with Tension Adjustable Bolt On Seat Upholstery)
 5" (12.7cm) 6" (15.2cm) 7" (17.8cm) 8" (20.3cm) 9" (22.9cm),.....**£306**
 ZERG2 Aero Z and ZRA Ergonomic Seat (Center of Gravity adjustment may be affected by this option. Not available with V-Front or Seat Taper, Seat Back that Locks When Folded, Desk Arms, Depth Adjustable Backrests, Fold Down Rigid Sideguards, Fendered Adjustable Sideguards, Solid Seat, Tension Adjustable Bolt On Seat Upholstery, or the Curve cushion.)
 6" (15.2cm) (Minimum Seat Depth B 12") 7" (17.8cm) (Minimum Seat Depth B 12")
 8" (20.3cm) (Minimum Seat Depth B 13")**£306**

FRONT SEAT HEIGHT C Measure from floor to top of seat tube at front edge of seat sling. ZRA adjusts up and down in .5" increments.

- ZFSH1 16" (40.6cm) 16.5" (41.9cm) 17" (43.2cm) 17.5" (44.5cm) 18" (45.7cm)
 18.5" (47cm) 19" (48.3cm) 19.5" (49.5cm) 20" (50.8cm) 20.5" (52.1cm)
 21" (53.3cm),.....**STD**





REAR SEAT HEIGHT D Measure from floor to top of seat tube at back post. ZRA adjusts up and down in 1/2" increments.

Rear Wheel Size	Rear Seat Height
20"	12.5" - 18.5" (31.8 - 47cm)
22"	14.0" - 19" (35.6 - 48.3cm)
24"	14.5" - 20" (36.8 - 50.8cm)
25"	14.5" - 20" (36.8 - 50.8cm)
26"	15.5" - 20" (39.4 - 50.8cm)

- ZRSH1 12.5" (31.8cm) 13" (33cm) 13.5" (34.3cm) 14" (35.6cm) 14.5" (36.8cm) 15" (38.1cm)
 15.5" (39.4cm) 16" (40.6cm) 16.5" (41.9cm) 17" (43.2cm) 17.5" (44.5cm) 18" (45.7cm) 18.5" (47cm)
 19" (48.3cm) 19.5" (49.5cm) 20" (50.8cm),.....**STD**

FRONT ANGLE L Measure from frame front to floor.

- ZFA1 70° 75° 80° 85° (Not available with 6" Front Wheels.) 90° (Not available with 5" or 6" Front Wheels.),.....**STD**
 ZFA2 Overall Frame Length M _____" Measure from front of backrest to Front Edge of Footrest.,.....**NCO**



International TiLite Z Series Order Form

TiFit: TAILORED INDIVIDUALIZED FRAME

SEAT TO FOOTEST E /FOOTESTS Measure from front edge of seat sling to top rear of footrest.

Standard Footrests Measurement must be at least 2" shorter than Front Seat Height C.
 13.5" (34.3cm) 14" (35.6cm) 14.5" (36.8cm) 15" (38.1cm) 15.5" (39.4cm) 16" (40.6cm)
 16.5" (41.9cm) 17" (43.2cm) 17.5" (44.5cm) 18" (45.7cm) 18.5" (47cm) 19" (48.3cm)

- ZFTR3 Angle Adjustable Footrest.....**STD**
- ZFTR1 Titanium Footrest w/Flat ABS Cover.....**£145**
- ZFTR2 Titanium Open Loop Footrest.....**£186**
- ZFTR12 Flip-Back Footrest (Not available with 10" or 11" Rear Seat Width A, Heavy Duty Frame, ABS Cover with 2" Sides or Front Roller.).....**£160**
- High Mount Footrest**
- 7"-9" 9"-11" 11"-13" Set at _____ (Use 1/2" increments)
- ZFTR4 Angle Adjustable High Mount Footrest.....**£75**

FOOTEST WIDTH I or V Select either a Standard Front End I if no Seat Taper T is selected or a V-Front End V if a Seat Taper T is selected.

STANDARD FRONT END I Find your Rear Seat Width A in the Chart below and select one of the options to the right. Measure from inside of front frame tube to inside of opposite front frame tube. (Not available with Front Seat Width (Seat Taper) T.)



- ZFE1 **NO V-Front Taper** (Required when both Heavy Duty Frame and Seat Taper T are selected.).....**STD**

REAR SEAT WIDTH A (Outside Measurement)	FOOTREST WIDTH I		
	No Taper <small>(2.5" Narrower than Rear Seat Width)</small>	1.5" Taper <small>(4" Narrower than Rear Seat Width)</small>	3.5" Taper <small>(6" Narrower than Rear Seat Width)</small>
10"	<input type="checkbox"/> 7.5"		
11"	<input type="checkbox"/> 8.5"		
12"	<input type="checkbox"/> 9.5"	<input type="checkbox"/> 8"	
13"	<input type="checkbox"/> 10.5"	<input type="checkbox"/> 9"	
14"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 10"	
15"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 11"	<input type="checkbox"/> 9"
16"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 10"
17"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 11"
18"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 12"
19"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 13"
20"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 14"

OR

V-FRONT END V Measure from inside of front frame tube to inside of opposite front frame tube 2-1/2" above footrest. (Required with Front Seat Taper T. Titanium Footrests have 1" of adjustment up and down. Not available with Ergonomic Seat frames.)

- ZFE2 V-Front Taper.....**£85**
- 9" (22.9cm) 10" (25.54cm) 11" (27.9cm)
- 12" (30.5cm) 13" (33cm) 14" (35.6cm)
- 15" (38.1cm)



International TiLite Z Series Order Form

TiFit: TAILORED INDIVIDUALIZED FRAME

BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

- ZBUP2 Velcro Adjustable (Not available with Fold-Down Push Handles.).....**STD**
Tension Adjustable by Straps (Required with Fold-Down Push Handles.)
 ZBUP3 Padded Nylon.....**£82**
 ZBUP5 AIR Breathable Upholstery (Not available with Ultrasuede Backrest Stripe.).....**£136**
 ZUPC3 Ultrasuede Backrest Stripe Black.....**£73**
 ZBUP0 Omit Back Upholstery.....**NCO**

SEAT BACK TYPE/PUSH HANDLES/SEAT BACK HEIGHT F Measure from top of back post to top of seat tube at rear of frame.

- ZSB1 Folding Aluminum Adjustable Height Back (Required for Heavy Duty Frame.).....**STD**
 ZSB2 Folding Titanium Adjustable Height Back and Rounded Rigidizer Bar (Titanium back will match Frame Finish. Not available with Bolt-On Push Handles.).....**£145**
 ZSB15 Non-Folding, Fixed Height Titanium with New Rounded Rigidizer Bar (Not Available with Aero Z, ZR or Transport Brackets. Titanium back will match Frame Finish.).....**£209**

TISHAFT™ BACK RELEASE BARS (Not available with Non-Folding, Fixed Height Titanium Back.)

- ZSB3 Aluminum.....**STD**
 ZSB6 New 3k Weave Carbon Fiber.....**£52**

PUSH HANDLES

- ZSB6 Integrated Push Handles.....**NCO**
 ZSB7 Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery.).....**£85**
 ZSB8 Bolt-On Push Handles (Includes 4" Deep Backrest Rigidizer Bar. Not available with Ti Back.).....**£145**
 ZSB9 Bolt-On Height-Adjustable Push Handles (Includes 4" Deep Backrest Rigidizer Bar. Not available with Ti Back.).....**£249**

FOLDING SEAT BACK HEIGHT F (ZSB1 and ZSB2)

Push Handles (Selected Above)	Short	Medium	Tall	Set at: (Use 1/2" increments)
None/ Bolt-On	<input type="checkbox"/> 9.5"-12.5" (24.1-31.7cm)	<input type="checkbox"/> 13"-17.5" (33.0-44.4cm)	<input type="checkbox"/> 17"-21.5" (33.0-44.4cm)	Set at: _____
Integrated	<input type="checkbox"/> 11.5"-13.5" (29.2-34.2cm)	<input type="checkbox"/> 14.5"-17.5" (36.8-44.4cm)	<input type="checkbox"/> 18.5"-21.5" (46.9-54.6cm)	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 14"-16" (35.5-40.6cm)	<input type="checkbox"/> 15.5"-19.5" (39.3-49.5cm)	<input type="checkbox"/> 19.5"-22" (49.5-55.8cm)	Set at: _____

NON-FOLDING, FIXED-HEIGHT F TITANIUM SEAT BACK with New Rounded Rigidizer Bar (ZSB15) (Titanium back will match Frame Finish. Required for fixed Side Guards. Not available with ZRA, Swing Away Tubular or Swing Away - Flip Back Tubular Armrests and Integrated Push Handles, Integrated Fold-Down Push Handles or Bolt-On Height Adjustable Push Handles.)

- 8" (20.3cm) 8.5" (21.6cm) 9" (22.9cm) 9.5" (24.1cm) 10" (25.4cm) 10.5" (26.7cm) 11" (27.9cm)
 11.5" (29.2cm) 12" (30.5cm) 12.5" (31.8cm) 13" (33cm) 13.5" (34.3cm) 14" (35.6) 14.5" (36.8cm)
 15" (38.1cm) 15.5" (39.4cm) 16" (40.6cm) 16.5" (41.9cm) 17" (43.2cm) 17.5" (44.5cm) 18" (45.7cm)
 18.5" (47cm) 19" (48.3cm) 19.5"(49.5cm) 20" (50.8cm)

AVAILABLE ON ALL BACKS

- ZSB11 Backrest that Locks when Folded (Not available with Ergonomic Seat.).....**STD**
 ZSB12 4" (10 cm) Deep Backrest Rigidizer Bar (Required with Heavy Duty Frame. Not available with a Non-Folding, Fixed Height Titanium Seat Back Type. 4" Deep Rigidizer Bar is NOT a rounded style rigidizer bar.).....**NCO**
 ZSB16 Depth Adjustable Back/Frame (Seat Depths B of 10"-18" allow 2" of rearward adjustment. Height Adjustable Push Handles have a limited height range with Depth Adjustable Back. Not Available with 19" or 20" Seat Depths B, Heavy Duty Frame, Tension Adjustable By Straps Seat Upholstery or Carbon Fiber Side Guards.).....**£152**

SEAT BACK ANGLE G Measure angle from front of back post to floor.

- ZSBA1 **Folding Back**
 80° 83° 86° 89° 92° 95° 98° 101°.....**STD**
 ZSBA1 **Non-Folding Back**
 80° 81° 82° 83° 84° 85° 86° 87° 88° 89° 90° 91° 92° 93° 94° 95° 96°
 97° 98° 99° 100° 101°.....**NCO**

CENTER OF GRAVITY H Measure from front of back post at the seat tube to the center of the rear axle. If no center of gravity is selected, chair will ship with the center of gravity in the rearward-most position.

- ZCOG1 0" .25" (.6cm) .5" (1.3cm) .75" (1.9cm) 1" (2.5cm) 1.25" (3.2cm) 1.5" (3.8cm) 1.75" (4.4cm)
 2" (5.1cm) 2.25" (5.7cm) 2.5" (6.4cm) 2.75" (7cm) 3" (7.6cm) 3.25" (8.3cm) 3.5" (8.9cm)
 3.75" (9.5cm) 4" (10.2cm) 4.25" (10.8cm) 4.5" (11.4cm) 4.75" (12.1cm) 5" (12.7cm).....**STD**

International TiLite Z Series Order Form

REAR WHEEL SPACING J /PLUGS Measure from outside of seat back post to inside of rear tire. 0.5" Rear Wheel Spacing on ZR not available with 0° or 2° Camber. ZRA 0.75" Rear Wheel Spacing not available with 0° or 2° Camber K. Ergonomic Seat selection may limit available center of gravity.

- ZRWS1 .75" (1.9cm) 1" (2.5cm) 1.25" (3.2cm) **STD**
- ZCBR12 Adjustable Rear Wheel Spacing Camber Plugs (Not available with Power Adaptable Frame.)..... **STD**
- ZCBR9 Additional Spacers to Increase Rear Wheel Spacing..... **NCO**
- ZCBR11 Zero Play Aluminum Camber Plugs (Not available with Power Adaptable Frame or more than 6° of Camber K.)..... **NCO**

TiFit: TAILORED INDIVIDUALIZED FRAME

SEAT UPHOLSTERY

- ZSUP1 Tension Adjustable Bolt On (Not available with Ergonomic Seat or Tapered Seat T. Black only.)..... **STD**
- ZSUP2 Tension Adjustable by Straps (Standard with Ergonomic Seat and Front Seat Width (Seat Taper) T of 1", 2" or 3". Black upholstery only.)..... **£127**
- ZSUP8 Aluminum Solid Seat (Not available with 10" or 11" Rear Seat Width A, Ergonomic Seat U requires CMD 016.)..... **£209**
- ZSUP0 Omit Seat..... **NCO**

SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

- ZCSH2 Comfort Curve Cushion (Stretch-Air cover) **£60**
- ZCSH11 ROHO Single Compartment (Fits seat widths 14" to 19") Profile: Low Mid High **£510**
- ZCSH14 ROHO Quadtro Select (Fits seat widths 14" to 19") Profile: Low Mid High **£510**
- ZCSH16 ROHO® Hybrid Select (Fits seat widths 14" to 18", and 20")..... **£635**

CUSTOM CONFIGURATION OPTIONS

CAMBER K /CAMBER TUBES

- ZCBR1 0° 2° 4° 6° (More than 0° of Camber K is not available with the Power Adaptable Frame. 8° of Camber K may be available through our Creative Mobility Division.)..... **STD**
- ZCBR15 Aluminum Camber Tube..... **STD**
- ZCBR3 Titanium Camber Tube (Not available with Power Adaptable Frame.)..... **£145**

FRONT WHEELS

- ZFW24 3"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire Silver Hub **OR** Black Hub..... **£130**
- ZFW14 4"x.75" Plastic Wheel w/Poly Tire (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)..... **NCO**
- ZFW3 4"x.75" Light-Up Micro-Caster (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)..... **£60**
- ZFW26 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.)..... **£70**
- ZFW27 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire Silver Hub **OR** Black Hub..... **£137**
- ZFW6 5"x1" Plastic Wheel w/Poly Tire (Not available with 90° Front Angle L, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)..... **NCO**
- ZFW28 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with 90° Front Angle L or Power Adaptable Frame.)..... **£70**
- ZFW30 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire Silver Hub **OR** Black Hub (Not available with 90° Front Angle L.)..... **£152**
- ZFW9 6"x1" Plastic Wheel w/Poly Tire (Not available with 85° or 90° Front Angle L, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)..... **NCO**
- ZFW32 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with 85° or 90° Front Angle L, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)..... **£75**
- ZFW33 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire Silver Hub **OR** Black Hub (Not available with 85° or 90° Front Angle L, Slipstream Fork, Glide™ Single-Sided Fork.)..... **£160**

PRECISION LOCK FRONT FORKS

- ZFK1 TiLite Standard Fork..... **STD**
- ZFK2 TiLite Slipstream Single-Sided Fork (Not available with Heavy Duty Frame or Power Adaptable Frame.)..... **£145**
- ZGL1 Out-Front Glide™ Suspension Fork (Black anodize only. Not available with Heavy Duty Frame.)..... **£358**
- Rider's Weight: 125 lbs. (56.69kg) 125-180 lbs. (56.69- 81.64kg) 180-265 lbs. (81.64- 120.2kg)
- ZTFS1 Titanium Fork Stem..... **£40**

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations.)

- ZRW14 TiLite SHADOW™ **STD**
- 20"(451)
- 22"(501)
- 24"(540)
- 25"(559) (Requires 15" minimum Seat Depth B.)
- 26"(590) (Requires 15" minimum Seat Depth B.)

PERFORMANCE REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations.)

- ZRW7 Spinergy SPOX Black Red Blue White..... **£431**
- 22"(501) - Std. Hub

International TiLite Z Series Order Form

- 24"(540) - Std. Hub
- 25"(559) - Std. Hub *(Requires 15" minimum Seat Depth B.)*
- 26"(590) - **Std. Hub** *(Requires 15" minimum Seat Depth B.)*

ZRW8 Spinergy LX Black Red Blue White Bright Orange Flamingo Pink Purple.....**£628**

- 22"(501) - Std. Hub
- 24"(540) - Std. Hub
- 25"(559) - Std. Hub *(Requires 15" minimum Seat Depth B.)*
- 26"(590) - **Std. Hub** *(Requires 15" minimum Seat Depth B.)*

ZRW9 **New Spinergy Design Lab LX Wheels** *(If same color is selected, anodized hubs and rims match color anodized package on page 1. Available with black spokes only.)*

- Gold Anodize Green Anodize Blue Anodize Red Anodize Purple Anodize Orange Anodize.....**£628**
- 22"(501) - Std. Hub
- 24"(540) - Std. Hub
- 25"(559) - Std. Hub *(Requires 15" minimum Seat Depth B.)*
- 26"(590) - Std. Hub *(Requires 15" minimum Seat Depth B.)*

ZRW16 **New Spinergy® Carbon Blade with 3k Weave Carbon Fiber** *(Not available with Heavy Duty Frame.)*.....**£953**

- 24"(540)
- 25"(559) *(Requires 15" minimum Seat Depth B.)*

ZRW16 Spinergy Carbon Blade *(Not available with Heavy Duty Frame, Natural-Fit or Natural-Fit LT with Standard Grip or Super Grip.)*.....**£953**

- 24"(540)
- 25"(559) *(Requires 15" minimum Seat Depth B.)*

ZRW17 X-Core Mag *(Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Sliver Anodized Handrims.)*.....**£82**

- 20"(451)
- 22"(501)
- 24"(540)

AXLES

- ZAXL1 Stainless Quick Release.....**STD**
- ZAXL2 Titanium Quick Release.....**£186**
- ZAXL3 Quad Release *(Not available with Mag or 26"(590) Spinergy Spox or LX Rear Wheels.)*.....**£115**

PNEUMATIC REAR TIRES

- ZRTR1 Treaded *(Not available with 25"(559) or 20" (451) Rear Wheels.)*.....**STD**
- ZRTR2 Airless Inserts *(Not available with Spinergy Carbon Blade Rear Wheels.)*.....**£70**
- ZRTR31 New MBL SpeedLite 20"(451), 22"(501), 24"(540), 25"(559) and 26"(590)**.....**STD**
- ZRTR30 24"(540), 25"(559) and 26"(590) Primo Sentinel.....**£97**
- ZRTR8 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution.....**£78**
- ZRTR9 24"(540) and 25"(559) Schwalbe® Slick.....**£78**
- ZRTR24 24"(540) Kenda Kobra Knobby.....**£60**

SOLID REAR TIRES *(Not available with Spinergy Carbon Blade Rear Wheels.)*

- ZRTR22 20"(451), 22"(501), 24"(540) and 25"(559) Primo Express.....**£82**
- ZRTR23 24"(540) and 25"(559) Primo Active.....**£100**

HANDRIMS

Handrims by OUT-FRONT®

- ZRIM1 Aluminum – Silver Anodized.....**STD**
- ZRIM21 Aluminum – Black Anodized *(Not available with 20"(451) Rear Wheels.)*.....**£40**
- ZRIM12 Titanium – Satin Finish *(Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.)*.....**£339**
- ZRIM10 The Surge (All-In-One Oval with Gription™ Strip, Larger Profile) *(Not available with 20"(451) Rear Wheels.)*.....**£272**
- ZRIM11 The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile) *(Not available with 20"(451) Rear Wheels.)*.....**£294**
- ZRIM20 Q-Grip (Exclusive High-Friction, Anti-Slip Coating) *(Not available with 20"(451) Rear Wheels.)*.....**£249**

The Natural-Fit® (Original, Larger Profile) *(Not available with 20"(451) Rear Wheels.)*

- ZRIM14 Standard Grip *(Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)*.....**£264**
- ZRIM15 Super Grip *(Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)*.....**£298**
- ZRIM16 Without Thumb Grip.....**£246**

The Natural-Fit® LT (Reduced Weight, Smaller Profile) *(Not available with 20"(451) Rear Wheels.)*

- ZRIM17 Standard Grip *(Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)*.....**£286**
- ZRIM18 Super Grip *(Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)*.....**£324**
- ZRIM19 Without Thumb Grip.....**£268**

Additional Handrim Options

- ZRIM5 Plastic Coated – Black.....**£97**
- ZRIM22 Omit Handrims.....**NCO**

Tab Length

- ZRIM7 Long Tabs.....**STD**
- ZRIM8 Short Tabs *(Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.)*.....**£15**

WHEEL LOCKS

International TiLite Z Series Order Form

Wheel Locks by OUT-FRONT®

ZWLK7	<input type="checkbox"/> Compact Composite Push to Lock (Available with .75" to 1.5" Rear Wheel Spacing) only.)	NCO
ZWLK8	<input type="checkbox"/> Compact Composite Pull to Lock (Available with .75" to 1.5" Rear Wheel Spacing) only.)	NCO
Z2WLK10	<input type="checkbox"/> Extension Handles for Compact Composite Push/Pull to Lock	£67

Additional Wheel Locks

ZWLK11	<input type="checkbox"/> New MBL Composite Scissor Lock	£108
ZWLK3	Uni-Lock Set At: <input type="checkbox"/> Push OR <input type="checkbox"/> Pull	£122
ZWLK4	<input type="checkbox"/> Extension Handles for Uni-Lock	£67
ZWLK15	<input type="checkbox"/> Grade Aids Push to Lock (Not available with solid rear tires)	£134

CUSTOM CONFIGURATION OPTIONS

SIDE GUARDS

Aluminum

ZSDG2	<input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std	£145
ZSDG3	<input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std	£152
ZSDG4	<input type="checkbox"/> Fold Down Rigid <input type="checkbox"/> Small OR <input type="checkbox"/> Std (Not available with Ergonomic Seat.)	£249
ZSDG6	<input type="checkbox"/> Fendered Adjustable (Available with 24" and 25" Rear Wheels. Requires Rear Seat Height D of 14.5" to 18.5".)	£249

New 3k Weave Carbon Fiber

ZSDG7	<input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std	£227
ZSDG8	<input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std	£234
ZSDG9	<input type="checkbox"/> Fold Down Rigid <input type="checkbox"/> Small OR <input type="checkbox"/> Std (Not available with Ergonomic Seat.)	£339
ZSDG11	<input type="checkbox"/> Fendered Adjustable (Available with 24" and 25" Rear Wheels. Requires Rear Seat Height D of 14.5" to 18.5".)	£335
ZSDG12	<input type="checkbox"/> Fendered Removable	£355

Available on Aluminum and Carbon Fiber Non-Fendered Side Guards

ZSDG13	<input type="checkbox"/> Side Guards Custom Cut to Profile of Rear Wheels	£231
--------	---	------

ARMRESTS (Desk Arms come with built-in Side Guards.)

ZARM1	Swing Away Tubular <input type="checkbox"/> 6.5" to 9.5" OR <input type="checkbox"/> 9.5" to 12.5" (Not Available with Heavy Duty Frame or Height Adjustable Push Handles.)	£160
ZARM4	<input type="checkbox"/> Swing Away Height Adjustable 8.5" to 11.5" - Flip Back Tubular (Not Available with Heavy Duty Frame.)	£231
ZARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard (Requires 14" minimum Seat Depth B. Not available with Ergonomic Seat.)	£249
	<input type="checkbox"/> 9" to 12" (22.9 to 30.5cm) OR <input type="checkbox"/> 12" to 15" (30.5-38.1cm)	
	<input type="checkbox"/> Standard Length Pad (10" (25.5cm)) OR <input type="checkbox"/> Full Length Pad (14" (35.5cm))	
ZARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 15" minimum Seat Depth B. Not Available with Ergonomic Seat.)	£272
	<input type="checkbox"/> 9" to 12" (22.9 to 30.5cm) OR <input type="checkbox"/> 12" to 15" (30.5-38.1cm)	
	<input type="checkbox"/> Standard Length Pad (10" (25.5cm)) OR <input type="checkbox"/> Full Length Pad (14" (35.5cm))	
Available Armrest Option		
ZARM5	<input type="checkbox"/> Covered in Black Ultrasuede (Not available with Desk Arms)	£82

SEAT BELTS

ZBLT1	<input type="checkbox"/> Velcro Adjustable, Black Nylon Webbing	£33
ZBLT2	<input type="checkbox"/> Bodypoint® 2" Push-Button Buckle (Auto-Style.)	£75
ZBLT4	<input type="checkbox"/> Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" or 20" Rear Seat Widths A.)	£137
ZBLT5	<input type="checkbox"/> Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" or 20" Rear Seat Widths A.)	£130
ZBLT8	<input type="checkbox"/> Q'Straint® WC-19 Wheelchair-Anchored Pelvic Belt (Recommended for use with Transit Tie Down. Not Available with ZR.)	£145
ZBLT10	<input type="checkbox"/> Bodypoint® Optional Band Clamp Mounting Hardware	£15
ZBLT16	<input type="checkbox"/> Bodypoint® EvoFlex® Pelvic Stabilizer Push Button	£142

CALF STRAP

ZCLF2	Velcro Adjustable	£26
	<input type="checkbox"/> Small (Recommended for 9" - 12" Footrest Widths.)	
	<input type="checkbox"/> Medium (Recommended for 11" - 15" Footrest Widths.)	
	<input type="checkbox"/> Large (Recommended for 14" - 17" Footrest Widths.)	
	<input type="checkbox"/> Xtra Large (Recommended for 16" - 19" Footrest Widths.)	
ZCLF3	Bodypoint® Padded Velcro Adjustable	£78
	<input type="checkbox"/> Small (Recommended for 8" - 11" Footrest Widths.)	
	<input type="checkbox"/> Medium (Recommended for 11" - 15" Footrest Widths.)	
	<input type="checkbox"/> Large (Recommended for 16" - 19" Footrest Widths.)	
ZCLF1	<input type="checkbox"/> Buckle Adjustable, Black Nylon Webbing	£18

International TiLite Z Series Order Form

Effective January 7, 2022 (New Items in Red)

MONO-TUBE RIGID FRAME

Page 8 of 8

BACKPACK

ZBCK3 Black Permobil Backpack.....£70

SEAT POUCH

ZPCH3 Black Permobil Seat Pouch.....£40

ZPCH4 Bodypoint® Mobility Bag.....£78

SPOKE GUARDS

ZGRD1 Clear Lexan with TiMotion Logo (Not available with Mag Rear Wheels.).....£97

CUSTOM CONFIGURATION OPTIONS

ANTI-TIPS

ZTIP1 Rear Aluminum - Black.....£119

ZTIP7 Black Aluminum Single.....£60

Mount on Left Side or Right Side.....NCO

ZTIP2 User Friendly Flip-up Pair.....£249

ZTIP8 User Friendly Flip-up Single.....£122

Mount on Left Side or Right Side.....NCO

ACCESSORIES

ZUEST1 Upper Extremity Support Tray 12" to 14" OR 13" to 15" Seat Widths A (Requires Armrests.).....£212

ZLUG1 Luggage Carriers.....£108

ZMPCT2 Neoprene Impact Guards with Ultrasuede Stripe Black.....£78

ZCH1 Crutch Holder.....£108

GENERAL

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Straint® is a registered trademark of Q'Straint Inc.

NOTES

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR